

Table of contents

\\Research Funded

RESEARCH

FHLIN

eegfmri_emoclips_s015_protocol

AAHead_Scout MBSIREPI_ref rest_MBSIREPI_acc_1.4s ep2d_bold_moco_for_coregistration AAHead_Scout MBSIREPI_ref emoclips1_MBSIREPI_acc_1.4s MBSIREPI_ref emoclips1_MBSIREPI_acc_1.4s ep2d_bold_moco_for_coregistration AAHead_Scout MBSIREPI_ref emoclips2_MBSIREPI_acc_1.4s ep2d_bold_moco_for_coregistration AAHead_Scout MBSIREPI_ref emoclips3_MBSIREPI_acc_1.4s ep2d_bold_moco_for_coregistration t1_mprage_sag_p2	iso1.0
--	--------

\\Research Funded\RESEARCH\FHLIN\eegefMRI_emoclips_s015_protocol\AAHead_Scout

TA: 0:14 PM: REF Voxel size: 1.6×1.6×1.6 mmPAT: 3 Rel. SNR: 1.00 : fl

Properties

Prio recon	On
Load images to viewer	On
Inline movie	Off
Auto store images	On
Load images to stamp segments	Off
Load images to graphic segments	On
Auto open inline display	Off
Auto close inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	Single measurement

Routine

Slab group	1
Slabs	1
Dist. factor	20 %
Position	L0.0 A10.0 H0.0 mm
Orientation	Sagittal
Phase enc. dir.	A >> P
Phase oversampling	0 %
Slice oversampling	0.0 %
Slices per slab	128
FoV read	260 mm
FoV phase	100.0 %
Slice thickness	1.6 mm
TR	3.15 ms
TE	1.37 ms
Averages	1
Concatenations	1
Filter	Prescan Normalize
Coil elements	HC1-7;NC1,2

Contrast - Common

TR	3.15 ms
TE	1.37 ms
Flip angle	8 deg

Contrast - Dynamic

Averages	1
Averaging mode	Short term
Reconstruction	Magnitude
Measurements	1

Resolution - Common

FoV read	260 mm
FoV phase	100.0 %
Slice thickness	1.6 mm
Base resolution	160
Phase resolution	100 %
Slice resolution	69 %
Phase partial Fourier	6/8
Slice partial Fourier	6/8
Trajectory	Cartesian

Resolution - iPAT

PAT mode	GRAPPA
Accel. factor PE	3
Ref. lines PE	24
Accel. factor 3D	1

Resolution - iPAT

Reference scan mode	Integrated
---------------------	------------

Resolution - Filter Image

Image Filter	Off
Distortion Corr.	Off
Prescan Normalize	On
Unfiltered images	Off
Normalize	Off
B1 filter	Off

Resolution - Filter Rawdata

Raw filter	Off
Elliptical filter	Off

Geometry - Common

Slab group	1
Slabs	1
Dist. factor	20 %
Position	L0.0 A10.0 H0.0 mm
Orientation	Sagittal
Phase enc. dir.	A >> P
Slice oversampling	0.0 %
Slices per slab	128
FoV read	260 mm
FoV phase	100.0 %
Slice thickness	1.6 mm
TR	3.15 ms
Multi-slice mode	Sequential
Series	Ascending
Concatenations	1

Geometry - AutoAlign

Slab group	1
Position	L0.0 A10.0 H0.0 mm
Orientation	Sagittal
Phase enc. dir.	A >> P
Initial Position	Isocenter
L	0.0 mm
P	0.0 mm
H	0.0 mm
Initial Rotation	0.00 deg
Initial Orientation	Transversal

Geometry - Tim Planning Suite

Set-n-Go Protocol	Off
Table position	H
Table position	0 mm
Inline Composing	Off

System - Miscellaneous

Positioning mode	REF
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Adaptive Combine
Save uncombined	Off
Matrix Optimization	Off

System - Miscellaneous

Coil Select Mode	Off - AutoCoilSelect
------------------	----------------------

System - Adjustments

B0 Shim mode	Tune up
B1 Shim mode	TrueForm
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto

System - Adjust Volume

Position	Isocenter
Orientation	Transversal
Rotation	0.00 deg
A >> P	263 mm
R >> L	350 mm
F >> H	350 mm
Reset	Off

System - pTx Volumes

B1 Shim mode	TrueForm
Excitation	Non-sel.

System - Tx/Rx

Frequency 1H	123.250023 MHz
Correction factor	1
Gain	Low
Img. Scale Cor.	1.000
Reset	Off
? Ref. amplitude 1H	0.000 V

Physio - PACE

Resp. control	Off
Concatenations	1

Inline - Common

Flip angle	8 deg
Measurements	1
Time to center	6.2 s

Inline - Inline

Subtract	Off
Measurements	1
StdDev	Off
Save original images	On

Inline - MIP

MIP-Sag	Off
MIP-Cor	Off
MIP-Tra	Off
MIP-Time	Off
Save original images	On

Inline - Composing

Inline Composing	Off
Distortion Corr.	Off

Inline - MapIt

Save original images	On
MapIt	None
Flip angle	8 deg
Measurements	1

Inline - MapIt

Contrasts	1
TR	3.15 ms
TE	1.37 ms

Sequence - Part 1

Introduction	On
Dimension	3D
Asymmetric echo	Weak
Contrasts	1
Multi-slice mode	Sequential
Bandwidth	540 Hz/Px

Sequence - Part 2

RF pulse type	Fast
Gradient mode	Normal
Excitation	Non-sel.
RF spoiling	On

Sequence - Assistant

Mode	Off
------	-----

\\Research Funded\RESEARCH\FHLIN\eeegfmri_emoclips_s015_protocol\MBSIREPI_ref

TA: 5.0 s PM: REF Voxel size: 5.0×5.0×4.0 mmRel. SNR: 1.00 : fl

Properties

Prio recon	Off
Load images to viewer	On
Inline movie	Off
Auto store images	On
Load images to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Auto close inline display	Off
Start measurement without further preparation	Off
Wait for user to start	On
Start measurements	Single measurement

Routine

Slice group	1
Slices	20
Dist. factor	25 %
Position	R5.9 A10.9 H6.9 mm
Orientation	T > C-13.7 > S-3.7
Phase enc. dir.	A >> P
AutoAlign	---
Phase oversampling	0 %
FoV read	210 mm
FoV phase	100.0 %
Slice thickness	4.0 mm
TR	50.00 ms
TE	27.60 ms
Averages	1
Concatenations	20
Filter	None
Coil elements	HC1-7

Contrast - Common

TR	50.00 ms
TE	27.60 ms
TD	0 ms
MTC	Off
Flip angle	30 deg
Fat suppr.	Fat sat.
Water suppr.	None

Contrast - Dynamic

Averages	1
Averaging mode	Short term
Reconstruction	Magnitude
Measurements	10
Pause after meas. 1	0.0 s
Pause after meas. 2	0.0 s
Pause after meas. 3	0.0 s
Pause after meas. 4	0.0 s
Pause after meas. 5	0.0 s
Pause after meas. 6	0.0 s
Pause after meas. 7	0.0 s
Pause after meas. 8	0.0 s
Pause after meas. 9	0.0 s
Multiple series	Off

Resolution - Common

FoV read	210 mm
FoV phase	100.0 %

Resolution - Common

Slice thickness	4.0 mm
Base resolution	42
Phase resolution	100 %
Phase partial Fourier	Off
Interpolation	Off

Resolution - Filter Image

Image Filter	Off
Distortion Corr.	Off
Prescan Normalize	Off
Normalize	Off
B1 filter	Off

Resolution - Filter Rawdata

Raw filter	Off
Elliptical filter	Off

Geometry - Common

Slice group	1
Slices	20
Dist. factor	25 %
Position	R5.9 A10.9 H6.9 mm
Orientation	T > C-13.7 > S-3.7
Phase enc. dir.	A >> P
FoV read	210 mm
FoV phase	100.0 %
Slice thickness	4.0 mm
TR	50.00 ms
Multi-slice mode	Sequential
Series	Ascending
Concatenations	20

Geometry - AutoAlign

Slice group	1
Position	R5.9 A10.9 H6.9 mm
Orientation	T > C-13.7 > S-3.7
Phase enc. dir.	A >> P
AutoAlign	---
Initial Position	R5.9 A10.9 H6.9
R	5.9 mm
A	10.9 mm
H	6.9 mm
Initial Rotation	-5.00 deg
Initial Orientation	T > C
T > C	-13.7
> S	-3.7

Geometry - Saturation

Fat suppr.	Fat sat.
Water suppr.	None
Special sat.	None

Geometry - Tim Planning Suite

Set-n-Go Protocol	Off
Table position	H
Table position	0 mm
Inline Composing	Off

System - Miscellaneous

Positioning mode	REF
------------------	-----

System - Miscellaneous

Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Adaptive Combine
Save uncombined	Off
Matrix Optimization	Off
AutoAlign	---
Coil Select Mode	On - Coil Memory

System - Adjustments

B0 Shim mode	Standard
B1 Shim mode	TrueForm
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto

System - Adjust Volume

Position	R5.9 A10.9 H6.9 mm
Orientation	T > C-13.7 > S-3.7
Rotation	-5.00 deg
A >> P	210 mm
R >> L	210 mm
F >> H	99 mm
Reset	Off

System - pTx Volumes

B1 Shim mode	TrueForm
--------------	----------

System - Tx/Rx

Frequency 1H	123.250023 MHz
Correction factor	1
Gain	High
Img. Scale Cor.	1.000
Reset	Off
? Ref. amplitude 1H	0.000 V

Physio - Signal1

1st Signal/Mode	None
TR	50.00 ms
Concatenations	20

Inline - Common

Subtract	Off
Measurements	10
StdDev	Off
Save original images	On

Inline - MIP

MIP-Sag	Off
MIP-Cor	Off
MIP-Tra	Off
MIP-Time	Off
Save original images	On

Inline - Composing

Inline Composing	Off
Distortion Corr.	Off

Sequence - Part 1

Introduction	Off
Dimension	2D
Contrasts	1
Multi-slice mode	Sequential
Bandwidth	3970 Hz/Px

Sequence - Part 2

Gradient mode	Fast
RF spoiling	Off

Sequence - Nuclei

TX/RX Nucleus	1H
TX/RX delta frequency	0 Hz
TX Nucleus	None
TX delta frequency	0 Hz
Coil elements	HC1-7

Sequence - Special

Reference Scan	On
Shift	3 FOV/*
Segment	2 #
crusher bfatSAT	0.0 ms
crusher afatSAT	0.4 ms
slewrates	22.00 T/m/s

Sequence - Assistant

Mode	Off
------	-----

\\Research Funded\RESEARCH\FHLIN\eeegfmri_emoclips_s015_protocol\rest_MBSIREPI_acc_1.4s

TA: 8:05 PM: FIX Voxel size: 5.0×5.0×4.0 mmRel. SNR: 1.00 : fl

Properties

Prio recon	Off
Load images to viewer	On
Inline movie	Off
Auto store images	On
Load images to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Auto close inline display	Off
Start measurement without further preparation	Off
Wait for user to start	On
Start measurements	Single measurement

Routine

Slice group	1
Slices	20
Dist. factor	25 %
Position	R5.9 A10.9 H6.9 mm
Orientation	T > C-13.7 > S-3.7
Phase enc. dir.	A >> P
AutoAlign	---
Phase oversampling	0 %
FoV read	210 mm
FoV phase	100.0 %
Slice thickness	4.0 mm
TR	50.00 ms
TE	27.60 ms
Averages	1
Concatenations	20
Filter	None
Coil elements	HC1-7

Contrast - Common

TR	50.00 ms
TE	27.60 ms
TD	0 ms
MTC	Off
Flip angle	30 deg
Fat suppr.	Fat sat.
Water suppr.	None

Contrast - Dynamic

Averages	1
Averaging mode	Short term
Reconstruction	Magnitude
Measurements	324
Pause after meas.	1.4 s
Multiple series	Off

Resolution - Common

FoV read	210 mm
FoV phase	100.0 %
Slice thickness	4.0 mm
Base resolution	42
Phase resolution	100 %
Phase partial Fourier	Off
Interpolation	Off

Resolution - Filter Image

Image Filter	Off
--------------	-----

Resolution - Filter Image

Distortion Corr.	Off
Prescan Normalize	Off
Normalize	Off
B1 filter	Off

Resolution - Filter Rawdata

Raw filter	Off
Elliptical filter	Off

Geometry - Common

Slice group	1
Slices	20
Dist. factor	25 %
Position	R5.9 A10.9 H6.9 mm
Orientation	T > C-13.7 > S-3.7
Phase enc. dir.	A >> P
FoV read	210 mm
FoV phase	100.0 %
Slice thickness	4.0 mm
TR	50.00 ms
Multi-slice mode	Sequential
Series	Ascending
Concatenations	20

Geometry - AutoAlign

Slice group	1
Position	R5.9 A10.9 H6.9 mm
Orientation	T > C-13.7 > S-3.7
Phase enc. dir.	A >> P
AutoAlign	---
Initial Position	R5.9 A10.9 H6.9
R	5.9 mm
A	10.9 mm
H	6.9 mm
Initial Rotation	-5.00 deg
Initial Orientation	T > C
T > C	-13.7
> S	-3.7

Geometry - Saturation

Fat suppr.	Fat sat.
Water suppr.	None
Special sat.	None

Geometry - Tim Planning Suite

Set-n-Go Protocol	Off
Table position	H
Table position	0 mm
Inline Composing	Off

System - Miscellaneous

Positioning mode	FIX
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Adaptive Combine
Save uncombined	Off

System - Miscellaneous

Matrix Optimization	Off
AutoAlign	---
Coil Select Mode	On - Coil Memory

System - Adjustments

B0 Shim mode	Standard
B1 Shim mode	TrueForm
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto

System - Adjust Volume

! Position	R5.9 A10.9 H6.9 mm
! Orientation	T > C-13.7 > S-3.7
! Rotation	-5.00 deg
! A >> P	210 mm
! R >> L	210 mm
! F >> H	99 mm
Reset	Off

System - pTx Volumes

B1 Shim mode	TrueForm
--------------	----------

System - Tx/Rx

Frequency 1H	123.250023 MHz
Correction factor	1
Gain	High
Img. Scale Cor.	1.000
Reset	Off
? Ref. amplitude 1H	0.000 V

Physio - Signal1

1st Signal/Mode	None
TR	50.00 ms
Concatenations	20

Inline - Common

Subtract	Off
Measurements	324
StdDev	Off
Save original images	On

Inline - MIP

MIP-Sag	Off
MIP-Cor	Off
MIP-Tra	Off
MIP-Time	Off
Save original images	On

Inline - Composing

Inline Composing	Off
Distortion Corr.	Off

Sequence - Part 1

Introduction	Off
Dimension	2D
Contrasts	1
Multi-slice mode	Sequential
Bandwidth	3970 Hz/Px

Sequence - Part 2

Gradient mode	Fast
RF spoiling	Off

Sequence - Nuclei

TX/RX Nucleus	1H
TX/RX delta frequency	0 Hz
TX Nucleus	None
TX delta frequency	0 Hz
Coil elements	HC1-7

Sequence - Special

Reference Scan	Off
Shift	3 FOV/*
Segment	2 #
crusher bfatSAT	0.0 ms
crusher afatSAT	0.4 ms
slewrate	22.00 T/m/s

Sequence - Assistant

Mode	Off
------	-----

\\Research Funded\RESEARCH\FHLIN\leegfmri_emoclips_s015_protocol\ep2d_bold_moco_for_coregistration

TA: 0:36 PM: FIX Voxel size: 3.3×3.3×4.0 mmPAT: Off Rel. SNR: 1.00 : epfid

Properties

Prio recon	Off
Load images to viewer	On
Inline movie	Off
Auto store images	On
Load images to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Auto close inline display	Off
Start measurement without further preparation	Off
Wait for user to start	On
Start measurements	Single measurement

Routine

Slice group	1
Slices	20
Dist. factor	25 %
Position	R5.9 A10.9 H6.9 mm
Orientation	T > C-13.7 > S-3.7
Phase enc. dir.	A >> P
AutoAlign	---
Phase oversampling	0 %
FoV read	210 mm
FoV phase	100.0 %
Slice thickness	4.0 mm
TR	2000 ms
TE	27.0 ms
Averages	1
Concatenations	1
Filter	None
Coil elements	HC1-7

Contrast - Common

TR	2000 ms
TE	27.0 ms
MTC	Off
Flip angle	60 deg
Fat suppr.	Fat sat.

Contrast - Dynamic

Averages	1
Averaging mode	Long term
Reconstruction	Magnitude
Measurements	15
Delay in TR	0 ms
Multiple series	Off

Resolution - Common

FoV read	210 mm
FoV phase	100.0 %
Slice thickness	4.0 mm
Base resolution	64
Phase resolution	100 %
Phase partial Fourier	Off
Interpolation	Off

Resolution - iPAT

Accel. mode	None
-------------	------

Resolution - Filter Image

Distortion Corr.	Off
Prescan Normalize	Off

Resolution - Filter Rawdata

Raw filter	Off
Elliptical filter	Off
Hamming	Off

Geometry - Common

Slice group	1
Slices	20
Dist. factor	25 %
Position	R5.9 A10.9 H6.9 mm
Orientation	T > C-13.7 > S-3.7
Phase enc. dir.	A >> P
FoV read	210 mm
FoV phase	100.0 %
Slice thickness	4.0 mm
TR	2000 ms
Multi-slice mode	Interleaved
Series	Interleaved
Concatenations	1

Geometry - AutoAlign

Slice group	1
Position	R5.9 A10.9 H6.9 mm
Orientation	T > C-13.7 > S-3.7
Phase enc. dir.	A >> P
AutoAlign	---
Initial Position	R5.9 A10.9 H6.9
R	5.9 mm
A	10.9 mm
H	6.9 mm
Initial Rotation	-5.00 deg
Initial Orientation	T > C
T > C	-13.7
> S	-3.7

Geometry - Saturation

Fat suppr.	Fat sat.
Special sat.	None

Geometry - Tim Planning Suite

Set-n-Go Protocol	Off
Table position	H
Table position	0 mm
Inline Composing	Off

System - Miscellaneous

Positioning mode	FIX
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Sum of Squares
Matrix Optimization	Performance
AutoAlign	---

System - Miscellaneous

Coil Select Mode	Off - AutoCoilSelect
------------------	----------------------

System - Adjustments

B0 Shim mode	Standard
B1 Shim mode	TrueForm
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto

System - Adjust Volume

! Position	R5.9 A10.9 H6.9 mm
! Orientation	T > C-13.7 > S-3.7
! Rotation	-5.00 deg
! A >> P	210 mm
! R >> L	210 mm
! F >> H	99 mm
Reset	Off

System - pTx Volumes

B1 Shim mode	TrueForm
Excitation	Standard

System - Tx/Rx

Frequency 1H	123.250023 MHz
Correction factor	1
Gain	High
Img. Scale Cor.	1.000
Reset	Off
? Ref. amplitude 1H	0.000 V

Physio - Signal1

1st Signal/Mode	None
TR	2000 ms
Concatenations	1

BOLD

GLM Statistics	On
Dynamic t-maps	Off
Ignore meas. at start	0
Ignore after transition	0
Model transition states	On
Temp. highpass filter	On
Threshold	4.00
Paradigm size	15
Meas[1]	Baseline
Meas[2]	Baseline
Meas[3]	Baseline
Meas[4]	Baseline
Meas[5]	Baseline
Meas[6]	Baseline
Meas[7]	Baseline
Meas[8]	Baseline
Meas[9]	Baseline
Meas[10]	Baseline
Meas[11]	Baseline
Meas[12]	Baseline
Meas[13]	Baseline
Meas[14]	Baseline
Meas[15]	Active
Motion correction	On
Spatial filter	Off
Measurements	15

BOLD

Delay in TR	0 ms
Multiple series	Off

Sequence - Part 1

Introduction	On
Multi-slice mode	Interleaved
Free echo spacing	Off
Echo spacing	0.49 ms
Bandwidth	2442 Hz/Px

Sequence - Part 2

EPI factor	64
RF pulse type	Normal
Gradient mode	Performance
Excitation	Standard

Sequence - pTX Pulses

\\Research Funded\RESEARCH\FHLIN\eegefMRI_emoclips_s015_protocol\AAHead_Scout

TA: 0:14 PM: REF Voxel size: 1.6×1.6×1.6 mmPAT: 3 Rel. SNR: 1.00 : fl

Properties

Prio recon	On
Load images to viewer	On
Inline movie	Off
Auto store images	On
Load images to stamp segments	Off
Load images to graphic segments	On
Auto open inline display	Off
Auto close inline display	Off
Start measurement without further preparation	Off
Wait for user to start	On
Start measurements	Single measurement

Routine

Slab group	1
Slabs	1
Dist. factor	20 %
Position	L0.0 A10.0 H0.0 mm
Orientation	Sagittal
Phase enc. dir.	A >> P
Phase oversampling	0 %
Slice oversampling	0.0 %
Slices per slab	128
FoV read	260 mm
FoV phase	100.0 %
Slice thickness	1.6 mm
TR	3.15 ms
TE	1.37 ms
Averages	1
Concatenations	1
Filter	Prescan Normalize
Coil elements	HC1-7;NC1,2

Contrast - Common

TR	3.15 ms
TE	1.37 ms
Flip angle	8 deg

Contrast - Dynamic

Averages	1
Averaging mode	Short term
Reconstruction	Magnitude
Measurements	1

Resolution - Common

FoV read	260 mm
FoV phase	100.0 %
Slice thickness	1.6 mm
Base resolution	160
Phase resolution	100 %
Slice resolution	69 %
Phase partial Fourier	6/8
Slice partial Fourier	6/8
Trajectory	Cartesian

Resolution - iPAT

PAT mode	GRAPPA
Accel. factor PE	3
Ref. lines PE	24
Accel. factor 3D	1

Resolution - iPAT

Reference scan mode	Integrated
---------------------	------------

Resolution - Filter Image

Image Filter	Off
Distortion Corr.	Off
Prescan Normalize	On
Unfiltered images	Off
Normalize	Off
B1 filter	Off

Resolution - Filter Rawdata

Raw filter	Off
Elliptical filter	Off

Geometry - Common

Slab group	1
Slabs	1
Dist. factor	20 %
Position	L0.0 A10.0 H0.0 mm
Orientation	Sagittal
Phase enc. dir.	A >> P
Slice oversampling	0.0 %
Slices per slab	128
FoV read	260 mm
FoV phase	100.0 %
Slice thickness	1.6 mm
TR	3.15 ms
Multi-slice mode	Sequential
Series	Ascending
Concatenations	1

Geometry - AutoAlign

Slab group	1
Position	L0.0 A10.0 H0.0 mm
Orientation	Sagittal
Phase enc. dir.	A >> P
Initial Position	Isocenter
L	0.0 mm
P	0.0 mm
H	0.0 mm
Initial Rotation	0.00 deg
Initial Orientation	Transversal

Geometry - Tim Planning Suite

Set-n-Go Protocol	Off
Table position	H
Table position	0 mm
Inline Composing	Off

System - Miscellaneous

Positioning mode	REF
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Adaptive Combine
Save uncombined	Off
Matrix Optimization	Off

System - Miscellaneous

Coil Select Mode	Off - AutoCoilSelect
------------------	----------------------

System - Adjustments

B0 Shim mode	Tune up
B1 Shim mode	TrueForm
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto

System - Adjust Volume

Position	Isocenter
Orientation	Transversal
Rotation	0.00 deg
A >> P	263 mm
R >> L	350 mm
F >> H	350 mm
Reset	Off

System - pTx Volumes

B1 Shim mode	TrueForm
Excitation	Non-sel.

System - Tx/Rx

Frequency 1H	123.250023 MHz
Correction factor	1
Gain	Low
Img. Scale Cor.	1.000
Reset	Off
? Ref. amplitude 1H	0.000 V

Physio - PACE

Resp. control	Off
Concatenations	1

Inline - Common

Flip angle	8 deg
Measurements	1
Time to center	6.2 s

Inline - Inline

Subtract	Off
Measurements	1
StdDev	Off
Save original images	On

Inline - MIP

MIP-Sag	Off
MIP-Cor	Off
MIP-Tra	Off
MIP-Time	Off
Save original images	On

Inline - Composing

Inline Composing	Off
Distortion Corr.	Off

Inline - MapIt

Save original images	On
MapIt	None
Flip angle	8 deg
Measurements	1

Inline - MapIt

Contrasts	1
TR	3.15 ms
TE	1.37 ms

Sequence - Part 1

Introduction	On
Dimension	3D
Asymmetric echo	Weak
Contrasts	1
Multi-slice mode	Sequential
Bandwidth	540 Hz/Px

Sequence - Part 2

RF pulse type	Fast
Gradient mode	Normal
Excitation	Non-sel.
RF spoiling	On

Sequence - Assistant

Mode	Off
------	-----

\\Research Funded\RESEARCH\FHLIN\eeegfmri_emoclips_s015_protocol\MBSIREPI_ref

TA: 5.0 s PM: FIX Voxel size: 5.0×5.0×4.0 mmRel. SNR: 1.00 : fl

Properties

Prio recon	Off
Load images to viewer	On
Inline movie	Off
Auto store images	On
Load images to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Auto close inline display	Off
Start measurement without further preparation	Off
Wait for user to start	On
Start measurements	Single measurement

Routine

Slice group	1
Slices	20
Dist. factor	25 %
Position	R5.9 A10.9 H6.9 mm
Orientation	T > C-13.7 > S-3.7
Phase enc. dir.	A >> P
AutoAlign	---
Phase oversampling	0 %
FoV read	210 mm
FoV phase	100.0 %
Slice thickness	4.0 mm
TR	50.00 ms
TE	27.60 ms
Averages	1
Concatenations	20
Filter	None
Coil elements	HC1-7

Contrast - Common

TR	50.00 ms
TE	27.60 ms
TD	0 ms
MTC	Off
Flip angle	30 deg
Fat suppr.	Fat sat.
Water suppr.	None

Contrast - Dynamic

Averages	1
Averaging mode	Short term
Reconstruction	Magnitude
Measurements	10
Pause after meas. 1	0.0 s
Pause after meas. 2	0.0 s
Pause after meas. 3	0.0 s
Pause after meas. 4	0.0 s
Pause after meas. 5	0.0 s
Pause after meas. 6	0.0 s
Pause after meas. 7	0.0 s
Pause after meas. 8	0.0 s
Pause after meas. 9	0.0 s
Multiple series	Off

Resolution - Common

FoV read	210 mm
FoV phase	100.0 %

Resolution - Common

Slice thickness	4.0 mm
Base resolution	42
Phase resolution	100 %
Phase partial Fourier	Off
Interpolation	Off

Resolution - Filter Image

Image Filter	Off
Distortion Corr.	Off
Prescan Normalize	Off
Normalize	Off
B1 filter	Off

Resolution - Filter Rawdata

Raw filter	Off
Elliptical filter	Off

Geometry - Common

Slice group	1
Slices	20
Dist. factor	25 %
Position	R5.9 A10.9 H6.9 mm
Orientation	T > C-13.7 > S-3.7
Phase enc. dir.	A >> P
FoV read	210 mm
FoV phase	100.0 %
Slice thickness	4.0 mm
TR	50.00 ms
Multi-slice mode	Sequential
Series	Ascending
Concatenations	20

Geometry - AutoAlign

Slice group	1
Position	R5.9 A10.9 H6.9 mm
Orientation	T > C-13.7 > S-3.7
Phase enc. dir.	A >> P
AutoAlign	---
Initial Position	R5.9 A10.9 H6.9
R	5.9 mm
A	10.9 mm
H	6.9 mm
Initial Rotation	-5.00 deg
Initial Orientation	T > C
T > C	-13.7
> S	-3.7

Geometry - Saturation

Fat suppr.	Fat sat.
Water suppr.	None
Special sat.	None

Geometry - Tim Planning Suite

Set-n-Go Protocol	Off
Table position	H
Table position	0 mm
Inline Composing	Off

System - Miscellaneous

Positioning mode	FIX
------------------	-----

System - Miscellaneous

Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Adaptive Combine
Save uncombined	Off
Matrix Optimization	Off
AutoAlign	---
Coil Select Mode	On - Coil Memory

System - Adjustments

B0 Shim mode	Standard
B1 Shim mode	TrueForm
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto

System - Adjust Volume

! Position	R5.9 A10.9 H6.9 mm
! Orientation	T > C-13.7 > S-3.7
! Rotation	-5.00 deg
! A >> P	210 mm
! R >> L	210 mm
! F >> H	99 mm
Reset	Off

System - pTx Volumes

B1 Shim mode	TrueForm
--------------	----------

System - Tx/Rx

Frequency 1H	123.250023 MHz
Correction factor	1
Gain	High
Img. Scale Cor.	1.000
Reset	Off
? Ref. amplitude 1H	0.000 V

Physio - Signal1

1st Signal/Mode	None
TR	50.00 ms
Concatenations	20

Inline - Common

Subtract	Off
Measurements	10
StdDev	Off
Save original images	On

Inline - MIP

MIP-Sag	Off
MIP-Cor	Off
MIP-Tra	Off
MIP-Time	Off
Save original images	On

Inline - Composing

Inline Composing	Off
Distortion Corr.	Off

Sequence - Part 1

Introduction	Off
Dimension	2D
Contrasts	1
Multi-slice mode	Sequential
Bandwidth	3970 Hz/Px

Sequence - Part 2

Gradient mode	Fast
RF spoiling	Off

Sequence - Nuclei

TX/RX Nucleus	1H
TX/RX delta frequency	0 Hz
TX Nucleus	None
TX delta frequency	0 Hz
Coil elements	HC1-7

Sequence - Special

Reference Scan	On
Shift	3 FOV/*
Segment	2 #
crusher bfatSAT	0.0 ms
crusher afatSAT	0.4 ms
slewrate	22.00 T/m/s

Sequence - Assistant

Mode	Off
------	-----

\\Research Funded\RESEARCH\FHLIN\eeegfmri_emoclips_s015_protocol\emoclips1_MBSIREPI_acc_1.4s

TA: 8:47 PM: FIX Voxel size: 5.0×5.0×4.0 mmRel. SNR: 1.00 : fl

Properties

Prio recon	Off
Load images to viewer	On
Inline movie	Off
Auto store images	On
Load images to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Auto close inline display	Off
Start measurement without further preparation	Off
Wait for user to start	On
Start measurements	Single measurement

Routine

Slice group	1
Slices	20
Dist. factor	25 %
Position	R5.9 A10.9 H6.9 mm
Orientation	T > C-13.7 > S-3.7
Phase enc. dir.	A >> P
AutoAlign	---
Phase oversampling	0 %
FoV read	210 mm
FoV phase	100.0 %
Slice thickness	4.0 mm
TR	50.00 ms
TE	27.60 ms
Averages	1
Concatenations	20
Filter	None
Coil elements	HC1-7

Contrast - Common

TR	50.00 ms
TE	27.60 ms
TD	0 ms
MTC	Off
Flip angle	30 deg
Fat suppr.	Fat sat.
Water suppr.	None

Contrast - Dynamic

Averages	1
Averaging mode	Short term
Reconstruction	Magnitude
Measurements	352
Pause after meas.	1.4 s
Multiple series	Off

Resolution - Common

FoV read	210 mm
FoV phase	100.0 %
Slice thickness	4.0 mm
Base resolution	42
Phase resolution	100 %
Phase partial Fourier	Off
Interpolation	Off

Resolution - Filter Image

Image Filter	Off
Distortion Corr.	Off
Prescan Normalize	Off
Normalize	Off
B1 filter	Off

Resolution - Filter Rawdata

Raw filter	Off
Elliptical filter	Off

Geometry - Common

Slice group	1
Slices	20
Dist. factor	25 %
Position	R5.9 A10.9 H6.9 mm
Orientation	T > C-13.7 > S-3.7
Phase enc. dir.	A >> P
FoV read	210 mm
FoV phase	100.0 %
Slice thickness	4.0 mm
TR	50.00 ms
Multi-slice mode	Sequential
Series	Ascending
Concatenations	20

Geometry - AutoAlign

Slice group	1
Position	R5.9 A10.9 H6.9 mm
Orientation	T > C-13.7 > S-3.7
Phase enc. dir.	A >> P
AutoAlign	---
Initial Position	R5.9 A10.9 H6.9
R	5.9 mm
A	10.9 mm
H	6.9 mm
Initial Rotation	-5.00 deg
Initial Orientation	T > C
T > C	-13.7
> S	-3.7

Geometry - Saturation

Fat suppr.	Fat sat.
Water suppr.	None
Special sat.	None

Geometry - Tim Planning Suite

Set-n-Go Protocol	Off
Table position	H
Table position	0 mm
Inline Composing	Off

System - Miscellaneous

Positioning mode	FIX
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H

System - Miscellaneous

Coil Combine Mode	Adaptive Combine
Save uncombined	Off
Matrix Optimization	Off
AutoAlign	---
Coil Select Mode	On - Coil Memory

System - Adjustments

B0 Shim mode	Standard
B1 Shim mode	TrueForm
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto

System - Adjust Volume

! Position	R5.9 A10.9 H6.9 mm
! Orientation	T > C-13.7 > S-3.7
! Rotation	-5.00 deg
! A >> P	210 mm
! R >> L	210 mm
! F >> H	99 mm
Reset	Off

System - pTx Volumes

B1 Shim mode	TrueForm
--------------	----------

System - Tx/Rx

Frequency 1H	123.250023 MHz
Correction factor	1
Gain	High
Img. Scale Cor.	1.000
Reset	Off
? Ref. amplitude 1H	0.000 V

Physio - Signal1

1st Signal/Mode	None
TR	50.00 ms
Concatenations	20

Inline - Common

Subtract	Off
Measurements	352
StdDev	Off
Save original images	On

Inline - MIP

MIP-Sag	Off
MIP-Cor	Off
MIP-Tra	Off
MIP-Time	Off
Save original images	On

Inline - Composing

Inline Composing	Off
Distortion Corr.	Off

Sequence - Part 1

Introduction	Off
Dimension	2D
Contrasts	1
Multi-slice mode	Sequential
Bandwidth	3970 Hz/Px

Sequence - Part 2

Gradient mode	Fast
RF spoiling	Off

Sequence - Nuclei

TX/RX Nucleus	1H
TX/RX delta frequency	0 Hz
TX Nucleus	None
TX delta frequency	0 Hz
Coil elements	HC1-7

Sequence - Special

Reference Scan	Off
Shift	3 FOV/*
Segment	2 #
crusher bfatSAT	0.0 ms
crusher afatSAT	0.4 ms
slewrate	22.00 T/m/s

Sequence - Assistant

Mode	Off
------	-----

\\Research Funded\RESEARCH\FHLIN\eeegfmri_emoclips_s015_protocol\MBSIREPI_ref

TA: 5.0 s PM: FIX Voxel size: 5.0×5.0×4.0 mmRel. SNR: 1.00 : fl

Properties

Prio recon	Off
Load images to viewer	On
Inline movie	Off
Auto store images	On
Load images to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Auto close inline display	Off
Start measurement without further preparation	Off
Wait for user to start	On
Start measurements	Single measurement

Routine

Slice group	1
Slices	20
Dist. factor	25 %
Position	R5.9 A10.9 H6.9 mm
Orientation	T > C-13.7 > S-3.7
Phase enc. dir.	A >> P
AutoAlign	---
Phase oversampling	0 %
FoV read	210 mm
FoV phase	100.0 %
Slice thickness	4.0 mm
TR	50.00 ms
TE	27.60 ms
Averages	1
Concatenations	20
Filter	None
Coil elements	HC1-7

Contrast - Common

TR	50.00 ms
TE	27.60 ms
TD	0 ms
MTC	Off
Flip angle	30 deg
Fat suppr.	Fat sat.
Water suppr.	None

Contrast - Dynamic

Averages	1
Averaging mode	Short term
Reconstruction	Magnitude
Measurements	10
Pause after meas. 1	0.0 s
Pause after meas. 2	0.0 s
Pause after meas. 3	0.0 s
Pause after meas. 4	0.0 s
Pause after meas. 5	0.0 s
Pause after meas. 6	0.0 s
Pause after meas. 7	0.0 s
Pause after meas. 8	0.0 s
Pause after meas. 9	0.0 s
Multiple series	Off

Resolution - Common

FoV read	210 mm
FoV phase	100.0 %

Resolution - Common

Slice thickness	4.0 mm
Base resolution	42
Phase resolution	100 %
Phase partial Fourier	Off
Interpolation	Off

Resolution - Filter Image

Image Filter	Off
Distortion Corr.	Off
Prescan Normalize	Off
Normalize	Off
B1 filter	Off

Resolution - Filter Rawdata

Raw filter	Off
Elliptical filter	Off

Geometry - Common

Slice group	1
Slices	20
Dist. factor	25 %
Position	R5.9 A10.9 H6.9 mm
Orientation	T > C-13.7 > S-3.7
Phase enc. dir.	A >> P
FoV read	210 mm
FoV phase	100.0 %
Slice thickness	4.0 mm
TR	50.00 ms
Multi-slice mode	Sequential
Series	Ascending
Concatenations	20

Geometry - AutoAlign

Slice group	1
Position	R5.9 A10.9 H6.9 mm
Orientation	T > C-13.7 > S-3.7
Phase enc. dir.	A >> P
AutoAlign	---
Initial Position	R5.9 A10.9 H6.9
R	5.9 mm
A	10.9 mm
H	6.9 mm
Initial Rotation	-5.00 deg
Initial Orientation	T > C
T > C	-13.7
> S	-3.7

Geometry - Saturation

Fat suppr.	Fat sat.
Water suppr.	None
Special sat.	None

Geometry - Tim Planning Suite

Set-n-Go Protocol	Off
Table position	H
Table position	0 mm
Inline Composing	Off

System - Miscellaneous

Positioning mode	FIX
------------------	-----

System - Miscellaneous

Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Adaptive Combine
Save uncombined	Off
Matrix Optimization	Off
AutoAlign	---
Coil Select Mode	On - Coil Memory

System - Adjustments

B0 Shim mode	Standard
B1 Shim mode	TrueForm
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto

System - Adjust Volume

! Position	R5.9 A10.9 H6.9 mm
! Orientation	T > C-13.7 > S-3.7
! Rotation	-5.00 deg
! A >> P	210 mm
! R >> L	210 mm
! F >> H	99 mm
Reset	Off

System - pTx Volumes

B1 Shim mode	TrueForm
--------------	----------

System - Tx/Rx

Frequency 1H	123.250023 MHz
Correction factor	1
Gain	High
Img. Scale Cor.	1.000
Reset	Off
? Ref. amplitude 1H	0.000 V

Physio - Signal1

1st Signal/Mode	None
TR	50.00 ms
Concatenations	20

Inline - Common

Subtract	Off
Measurements	10
StdDev	Off
Save original images	On

Inline - MIP

MIP-Sag	Off
MIP-Cor	Off
MIP-Tra	Off
MIP-Time	Off
Save original images	On

Inline - Composing

Inline Composing	Off
Distortion Corr.	Off

Sequence - Part 1

Introduction	Off
Dimension	2D
Contrasts	1
Multi-slice mode	Sequential
Bandwidth	3970 Hz/Px

Sequence - Part 2

Gradient mode	Fast
RF spoiling	Off

Sequence - Nuclei

TX/RX Nucleus	1H
TX/RX delta frequency	0 Hz
TX Nucleus	None
TX delta frequency	0 Hz
Coil elements	HC1-7

Sequence - Special

Reference Scan	On
Shift	3 FOV/*
Segment	2 #
crusher bfatSAT	0.0 ms
crusher afatSAT	0.4 ms
slewrates	22.00 T/m/s

Sequence - Assistant

Mode	Off
------	-----

\\Research Funded\RESEARCH\FHLIN\eeegfmri_emoclips_s015_protocol\emoclips1_MBSIREPI_acc_1.4s

TA: 8:47 PM: FIX Voxel size: 5.0×5.0×4.0 mmRel. SNR: 1.00 : fl

Properties

Prio recon	Off
Load images to viewer	On
Inline movie	Off
Auto store images	On
Load images to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Auto close inline display	Off
Start measurement without further preparation	Off
Wait for user to start	On
Start measurements	Single measurement

Routine

Slice group	1
Slices	20
Dist. factor	25 %
Position	R5.9 A10.9 H6.9 mm
Orientation	T > C-13.7 > S-3.7
Phase enc. dir.	A >> P
AutoAlign	---
Phase oversampling	0 %
FoV read	210 mm
FoV phase	100.0 %
Slice thickness	4.0 mm
TR	50.00 ms
TE	27.60 ms
Averages	1
Concatenations	20
Filter	None
Coil elements	HC1-7

Contrast - Common

TR	50.00 ms
TE	27.60 ms
TD	0 ms
MTC	Off
Flip angle	30 deg
Fat suppr.	Fat sat.
Water suppr.	None

Contrast - Dynamic

Averages	1
Averaging mode	Short term
Reconstruction	Magnitude
Measurements	352
Pause after meas.	1.4 s
Multiple series	Off

Resolution - Common

FoV read	210 mm
FoV phase	100.0 %
Slice thickness	4.0 mm
Base resolution	42
Phase resolution	100 %
Phase partial Fourier	Off
Interpolation	Off

Resolution - Filter Image

Image Filter	Off
Distortion Corr.	Off
Prescan Normalize	Off
Normalize	Off
B1 filter	Off

Resolution - Filter Rawdata

Raw filter	Off
Elliptical filter	Off

Geometry - Common

Slice group	1
Slices	20
Dist. factor	25 %
Position	R5.9 A10.9 H6.9 mm
Orientation	T > C-13.7 > S-3.7
Phase enc. dir.	A >> P
FoV read	210 mm
FoV phase	100.0 %
Slice thickness	4.0 mm
TR	50.00 ms
Multi-slice mode	Sequential
Series	Ascending
Concatenations	20

Geometry - AutoAlign

Slice group	1
Position	R5.9 A10.9 H6.9 mm
Orientation	T > C-13.7 > S-3.7
Phase enc. dir.	A >> P
AutoAlign	---
Initial Position	R5.9 A10.9 H6.9
R	5.9 mm
A	10.9 mm
H	6.9 mm
Initial Rotation	-5.00 deg
Initial Orientation	T > C
T > C	-13.7
> S	-3.7

Geometry - Saturation

Fat suppr.	Fat sat.
Water suppr.	None
Special sat.	None

Geometry - Tim Planning Suite

Set-n-Go Protocol	Off
Table position	H
Table position	0 mm
Inline Composing	Off

System - Miscellaneous

Positioning mode	FIX
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H

System - Miscellaneous

Coil Combine Mode	Adaptive Combine
Save uncombined	Off
Matrix Optimization	Off
AutoAlign	---
Coil Select Mode	On - Coil Memory

System - Adjustments

B0 Shim mode	Standard
B1 Shim mode	TrueForm
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto

System - Adjust Volume

! Position	R5.9 A10.9 H6.9 mm
! Orientation	T > C-13.7 > S-3.7
! Rotation	-5.00 deg
! A >> P	210 mm
! R >> L	210 mm
! F >> H	99 mm
Reset	Off

System - pTx Volumes

B1 Shim mode	TrueForm
--------------	----------

System - Tx/Rx

Frequency 1H	123.250023 MHz
Correction factor	1
Gain	High
Img. Scale Cor.	1.000
Reset	Off
? Ref. amplitude 1H	0.000 V

Physio - Signal1

1st Signal/Mode	None
TR	50.00 ms
Concatenations	20

Inline - Common

Subtract	Off
Measurements	352
StdDev	Off
Save original images	On

Inline - MIP

MIP-Sag	Off
MIP-Cor	Off
MIP-Tra	Off
MIP-Time	Off
Save original images	On

Inline - Composing

Inline Composing	Off
Distortion Corr.	Off

Sequence - Part 1

Introduction	Off
Dimension	2D
Contrasts	1
Multi-slice mode	Sequential
Bandwidth	3970 Hz/Px

Sequence - Part 2

Gradient mode	Fast
RF spoiling	Off

Sequence - Nuclei

TX/RX Nucleus	1H
TX/RX delta frequency	0 Hz
TX Nucleus	None
TX delta frequency	0 Hz
Coil elements	HC1-7

Sequence - Special

Reference Scan	Off
Shift	3 FOV/*
Segment	2 #
crusher bfatSAT	0.0 ms
crusher afatSAT	0.4 ms
slewrate	22.00 T/m/s

Sequence - Assistant

Mode	Off
------	-----

\\Research Funded\RESEARCH\FHLIN\leegfmri_emoclips_s015_protocol\ep2d_bold_moco_for_coregistration

TA: 0:36 PM: FIX Voxel size: 3.3×3.3×4.0 mmPAT: Off Rel. SNR: 1.00 : epfid

Properties

Prio recon	Off
Load images to viewer	On
Inline movie	Off
Auto store images	On
Load images to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Auto close inline display	Off
Start measurement without further preparation	Off
Wait for user to start	On
Start measurements	Single measurement

Routine

Slice group	1
Slices	20
Dist. factor	25 %
Position	R5.9 A10.9 H6.9 mm
Orientation	T > C-13.7 > S-3.7
Phase enc. dir.	A >> P
AutoAlign	---
Phase oversampling	0 %
FoV read	210 mm
FoV phase	100.0 %
Slice thickness	4.0 mm
TR	2000 ms
TE	27.0 ms
Averages	1
Concatenations	1
Filter	None
Coil elements	HC1-7

Contrast - Common

TR	2000 ms
TE	27.0 ms
MTC	Off
Flip angle	60 deg
Fat suppr.	Fat sat.

Contrast - Dynamic

Averages	1
Averaging mode	Long term
Reconstruction	Magnitude
Measurements	15
Delay in TR	0 ms
Multiple series	Off

Resolution - Common

FoV read	210 mm
FoV phase	100.0 %
Slice thickness	4.0 mm
Base resolution	64
Phase resolution	100 %
Phase partial Fourier	Off
Interpolation	Off

Resolution - iPAT

Accel. mode	None
-------------	------

Resolution - Filter Image

Distortion Corr.	Off
Prescan Normalize	Off

Resolution - Filter Rawdata

Raw filter	Off
Elliptical filter	Off
Hamming	Off

Geometry - Common

Slice group	1
Slices	20
Dist. factor	25 %
Position	R5.9 A10.9 H6.9 mm
Orientation	T > C-13.7 > S-3.7
Phase enc. dir.	A >> P
FoV read	210 mm
FoV phase	100.0 %
Slice thickness	4.0 mm
TR	2000 ms
Multi-slice mode	Interleaved
Series	Interleaved
Concatenations	1

Geometry - AutoAlign

Slice group	1
Position	R5.9 A10.9 H6.9 mm
Orientation	T > C-13.7 > S-3.7
Phase enc. dir.	A >> P
AutoAlign	---
Initial Position	R5.9 A10.9 H6.9
R	5.9 mm
A	10.9 mm
H	6.9 mm
Initial Rotation	-5.00 deg
Initial Orientation	T > C
T > C	-13.7
> S	-3.7

Geometry - Saturation

Fat suppr.	Fat sat.
Special sat.	None

Geometry - Tim Planning Suite

Set-n-Go Protocol	Off
Table position	H
Table position	0 mm
Inline Composing	Off

System - Miscellaneous

Positioning mode	FIX
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Sum of Squares
Matrix Optimization	Performance
AutoAlign	---

System - Miscellaneous

Coil Select Mode	Off - AutoCoilSelect
------------------	----------------------

System - Adjustments

B0 Shim mode	Standard
B1 Shim mode	TrueForm
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto

System - Adjust Volume

! Position	R5.9 A10.9 H6.9 mm
! Orientation	T > C-13.7 > S-3.7
! Rotation	-5.00 deg
! A >> P	210 mm
! R >> L	210 mm
! F >> H	99 mm
Reset	Off

System - pTx Volumes

B1 Shim mode	TrueForm
Excitation	Standard

System - Tx/Rx

Frequency 1H	123.250023 MHz
Correction factor	1
Gain	High
Img. Scale Cor.	1.000
Reset	Off
? Ref. amplitude 1H	0.000 V

Physio - Signal1

1st Signal/Mode	None
TR	2000 ms
Concatenations	1

BOLD

GLM Statistics	On
Dynamic t-maps	Off
Ignore meas. at start	0
Ignore after transition	0
Model transition states	On
Temp. highpass filter	On
Threshold	4.00
Paradigm size	15
Meas[1]	Baseline
Meas[2]	Baseline
Meas[3]	Baseline
Meas[4]	Baseline
Meas[5]	Baseline
Meas[6]	Baseline
Meas[7]	Baseline
Meas[8]	Baseline
Meas[9]	Baseline
Meas[10]	Baseline
Meas[11]	Baseline
Meas[12]	Baseline
Meas[13]	Baseline
Meas[14]	Baseline
Meas[15]	Active
Motion correction	On
Spatial filter	Off
Measurements	15

BOLD

Delay in TR	0 ms
Multiple series	Off

Sequence - Part 1

Introduction	On
Multi-slice mode	Interleaved
Free echo spacing	Off
Echo spacing	0.49 ms
Bandwidth	2442 Hz/Px

Sequence - Part 2

EPI factor	64
RF pulse type	Normal
Gradient mode	Performance
Excitation	Standard

Sequence - pTX Pulses

\\Research Funded\RESEARCH\FHLIN\eegefMRI_emoclips_s015_protocol\AAHead_Scout

TA: 0:14 PM: REF Voxel size: 1.6×1.6×1.6 mmPAT: 3 Rel. SNR: 1.00 : fl

Properties

Prio recon	On
Load images to viewer	On
Inline movie	Off
Auto store images	On
Load images to stamp segments	Off
Load images to graphic segments	On
Auto open inline display	Off
Auto close inline display	Off
Start measurement without further preparation	Off
Wait for user to start	On
Start measurements	Single measurement

Routine

Slab group	1
Slabs	1
Dist. factor	20 %
Position	L0.0 A10.0 H0.0 mm
Orientation	Sagittal
Phase enc. dir.	A >> P
Phase oversampling	0 %
Slice oversampling	0.0 %
Slices per slab	128
FoV read	260 mm
FoV phase	100.0 %
Slice thickness	1.6 mm
TR	3.15 ms
TE	1.37 ms
Averages	1
Concatenations	1
Filter	Prescan Normalize
Coil elements	HC1-7;NC1,2

Contrast - Common

TR	3.15 ms
TE	1.37 ms
Flip angle	8 deg

Contrast - Dynamic

Averages	1
Averaging mode	Short term
Reconstruction	Magnitude
Measurements	1

Resolution - Common

FoV read	260 mm
FoV phase	100.0 %
Slice thickness	1.6 mm
Base resolution	160
Phase resolution	100 %
Slice resolution	69 %
Phase partial Fourier	6/8
Slice partial Fourier	6/8
Trajectory	Cartesian

Resolution - iPAT

PAT mode	GRAPPA
Accel. factor PE	3
Ref. lines PE	24
Accel. factor 3D	1

Resolution - iPAT

Reference scan mode	Integrated
---------------------	------------

Resolution - Filter Image

Image Filter	Off
Distortion Corr.	Off
Prescan Normalize	On
Unfiltered images	Off
Normalize	Off
B1 filter	Off

Resolution - Filter Rawdata

Raw filter	Off
Elliptical filter	Off

Geometry - Common

Slab group	1
Slabs	1
Dist. factor	20 %
Position	L0.0 A10.0 H0.0 mm
Orientation	Sagittal
Phase enc. dir.	A >> P
Slice oversampling	0.0 %
Slices per slab	128
FoV read	260 mm
FoV phase	100.0 %
Slice thickness	1.6 mm
TR	3.15 ms
Multi-slice mode	Sequential
Series	Ascending
Concatenations	1

Geometry - AutoAlign

Slab group	1
Position	L0.0 A10.0 H0.0 mm
Orientation	Sagittal
Phase enc. dir.	A >> P
Initial Position	Isocenter
L	0.0 mm
P	0.0 mm
H	0.0 mm
Initial Rotation	0.00 deg
Initial Orientation	Transversal

Geometry - Tim Planning Suite

Set-n-Go Protocol	Off
Table position	H
Table position	0 mm
Inline Composing	Off

System - Miscellaneous

Positioning mode	REF
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Adaptive Combine
Save uncombined	Off
Matrix Optimization	Off

System - Miscellaneous

Coil Select Mode	Off - AutoCoilSelect
------------------	----------------------

System - Adjustments

B0 Shim mode	Tune up
B1 Shim mode	TrueForm
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto

System - Adjust Volume

Position	Isocenter
Orientation	Transversal
Rotation	0.00 deg
A >> P	263 mm
R >> L	350 mm
F >> H	350 mm
Reset	Off

System - pTx Volumes

B1 Shim mode	TrueForm
Excitation	Non-sel.

System - Tx/Rx

Frequency 1H	123.250023 MHz
Correction factor	1
Gain	Low
Img. Scale Cor.	1.000
Reset	Off
? Ref. amplitude 1H	0.000 V

Physio - PACE

Resp. control	Off
Concatenations	1

Inline - Common

Flip angle	8 deg
Measurements	1
Time to center	6.2 s

Inline - Inline

Subtract	Off
Measurements	1
StdDev	Off
Save original images	On

Inline - MIP

MIP-Sag	Off
MIP-Cor	Off
MIP-Tra	Off
MIP-Time	Off
Save original images	On

Inline - Composing

Inline Composing	Off
Distortion Corr.	Off

Inline - MapIt

Save original images	On
MapIt	None
Flip angle	8 deg
Measurements	1

Inline - MapIt

Contrasts	1
TR	3.15 ms
TE	1.37 ms

Sequence - Part 1

Introduction	On
Dimension	3D
Asymmetric echo	Weak
Contrasts	1
Multi-slice mode	Sequential
Bandwidth	540 Hz/Px

Sequence - Part 2

RF pulse type	Fast
Gradient mode	Normal
Excitation	Non-sel.
RF spoiling	On

Sequence - Assistant

Mode	Off
------	-----

\\Research Funded\RESEARCH\FHLIN\eeegfmri_emoclips_s015_protocol\MBSIREPI_ref

TA: 5.0 s PM: FIX Voxel size: 5.0×5.0×4.0 mmRel. SNR: 1.00 : fl

Properties

Prio recon	Off
Load images to viewer	On
Inline movie	Off
Auto store images	On
Load images to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Auto close inline display	Off
Start measurement without further preparation	Off
Wait for user to start	On
Start measurements	Single measurement

Routine

Slice group	1
Slices	20
Dist. factor	25 %
Position	R5.9 A10.9 H6.9 mm
Orientation	T > C-13.7 > S-3.7
Phase enc. dir.	A >> P
AutoAlign	---
Phase oversampling	0 %
FoV read	210 mm
FoV phase	100.0 %
Slice thickness	4.0 mm
TR	50.00 ms
TE	27.60 ms
Averages	1
Concatenations	20
Filter	None
Coil elements	HC1-7

Contrast - Common

TR	50.00 ms
TE	27.60 ms
TD	0 ms
MTC	Off
Flip angle	30 deg
Fat suppr.	Fat sat.
Water suppr.	None

Contrast - Dynamic

Averages	1
Averaging mode	Short term
Reconstruction	Magnitude
Measurements	10
Pause after meas. 1	0.0 s
Pause after meas. 2	0.0 s
Pause after meas. 3	0.0 s
Pause after meas. 4	0.0 s
Pause after meas. 5	0.0 s
Pause after meas. 6	0.0 s
Pause after meas. 7	0.0 s
Pause after meas. 8	0.0 s
Pause after meas. 9	0.0 s
Multiple series	Off

Resolution - Common

FoV read	210 mm
FoV phase	100.0 %

Resolution - Common

Slice thickness	4.0 mm
Base resolution	42
Phase resolution	100 %
Phase partial Fourier	Off
Interpolation	Off

Resolution - Filter Image

Image Filter	Off
Distortion Corr.	Off
Prescan Normalize	Off
Normalize	Off
B1 filter	Off

Resolution - Filter Rawdata

Raw filter	Off
Elliptical filter	Off

Geometry - Common

Slice group	1
Slices	20
Dist. factor	25 %
Position	R5.9 A10.9 H6.9 mm
Orientation	T > C-13.7 > S-3.7
Phase enc. dir.	A >> P
FoV read	210 mm
FoV phase	100.0 %
Slice thickness	4.0 mm
TR	50.00 ms
Multi-slice mode	Sequential
Series	Ascending
Concatenations	20

Geometry - AutoAlign

Slice group	1
Position	R5.9 A10.9 H6.9 mm
Orientation	T > C-13.7 > S-3.7
Phase enc. dir.	A >> P
AutoAlign	---
Initial Position	R5.9 A10.9 H6.9
R	5.9 mm
A	10.9 mm
H	6.9 mm
Initial Rotation	-5.00 deg
Initial Orientation	T > C
T > C	-13.7
> S	-3.7

Geometry - Saturation

Fat suppr.	Fat sat.
Water suppr.	None
Special sat.	None

Geometry - Tim Planning Suite

Set-n-Go Protocol	Off
Table position	H
Table position	0 mm
Inline Composing	Off

System - Miscellaneous

Positioning mode	FIX
------------------	-----

System - Miscellaneous

Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Adaptive Combine
Save uncombined	Off
Matrix Optimization	Off
AutoAlign	---
Coil Select Mode	On - Coil Memory

System - Adjustments

B0 Shim mode	Standard
B1 Shim mode	TrueForm
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto

System - Adjust Volume

! Position	R5.9 A10.9 H6.9 mm
! Orientation	T > C-13.7 > S-3.7
! Rotation	-5.00 deg
! A >> P	210 mm
! R >> L	210 mm
! F >> H	99 mm
Reset	Off

System - pTx Volumes

B1 Shim mode	TrueForm
--------------	----------

System - Tx/Rx

Frequency 1H	123.250023 MHz
Correction factor	1
Gain	High
Img. Scale Cor.	1.000
Reset	Off
? Ref. amplitude 1H	0.000 V

Physio - Signal1

1st Signal/Mode	None
TR	50.00 ms
Concatenations	20

Inline - Common

Subtract	Off
Measurements	10
StdDev	Off
Save original images	On

Inline - MIP

MIP-Sag	Off
MIP-Cor	Off
MIP-Tra	Off
MIP-Time	Off
Save original images	On

Inline - Composing

Inline Composing	Off
Distortion Corr.	Off

Sequence - Part 1

Introduction	Off
Dimension	2D
Contrasts	1
Multi-slice mode	Sequential
Bandwidth	3970 Hz/Px

Sequence - Part 2

Gradient mode	Fast
RF spoiling	Off

Sequence - Nuclei

TX/RX Nucleus	1H
TX/RX delta frequency	0 Hz
TX Nucleus	None
TX delta frequency	0 Hz
Coil elements	HC1-7

Sequence - Special

Reference Scan	On
Shift	3 FOV/*
Segment	2 #
crusher bfatSAT	0.0 ms
crusher afatSAT	0.4 ms
slewrate	22.00 T/m/s

Sequence - Assistant

Mode	Off
------	-----

\\Research Funded\RESEARCH\FHLIN\eeegfmri_emoclips_s015_protocol\emoclips2_MBSIREPI_acc_1.4s

TA: 8:06 PM: FIX Voxel size: 5.0×5.0×4.0 mmRel. SNR: 1.00 : fl

Properties

Prio recon	Off
Load images to viewer	On
Inline movie	Off
Auto store images	On
Load images to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Auto close inline display	Off
Start measurement without further preparation	Off
Wait for user to start	On
Start measurements	Single measurement

Routine

Slice group	1
Slices	20
Dist. factor	25 %
Position	R5.9 A10.9 H6.9 mm
Orientation	T > C-13.7 > S-3.7
Phase enc. dir.	A >> P
AutoAlign	---
Phase oversampling	0 %
FoV read	210 mm
FoV phase	100.0 %
Slice thickness	4.0 mm
TR	50.00 ms
TE	27.60 ms
Averages	1
Concatenations	20
Filter	None
Coil elements	HC1-7

Contrast - Common

TR	50.00 ms
TE	27.60 ms
TD	0 ms
MTC	Off
Flip angle	30 deg
Fat suppr.	Fat sat.
Water suppr.	None

Contrast - Dynamic

Averages	1
Averaging mode	Short term
Reconstruction	Magnitude
Measurements	325
Pause after meas.	1.4 s
Multiple series	Off

Resolution - Common

FoV read	210 mm
FoV phase	100.0 %
Slice thickness	4.0 mm
Base resolution	42
Phase resolution	100 %
Phase partial Fourier	Off
Interpolation	Off

Resolution - Filter Image

Image Filter	Off
Distortion Corr.	Off
Prescan Normalize	Off
Normalize	Off
B1 filter	Off

Resolution - Filter Rawdata

Raw filter	Off
Elliptical filter	Off

Geometry - Common

Slice group	1
Slices	20
Dist. factor	25 %
Position	R5.9 A10.9 H6.9 mm
Orientation	T > C-13.7 > S-3.7
Phase enc. dir.	A >> P
FoV read	210 mm
FoV phase	100.0 %
Slice thickness	4.0 mm
TR	50.00 ms
Multi-slice mode	Sequential
Series	Ascending
Concatenations	20

Geometry - AutoAlign

Slice group	1
Position	R5.9 A10.9 H6.9 mm
Orientation	T > C-13.7 > S-3.7
Phase enc. dir.	A >> P
AutoAlign	---
Initial Position	R5.9 A10.9 H6.9
R	5.9 mm
A	10.9 mm
H	6.9 mm
Initial Rotation	-5.00 deg
Initial Orientation	T > C
T > C	-13.7
> S	-3.7

Geometry - Saturation

Fat suppr.	Fat sat.
Water suppr.	None
Special sat.	None

Geometry - Tim Planning Suite

Set-n-Go Protocol	Off
Table position	H
Table position	0 mm
Inline Composing	Off

System - Miscellaneous

Positioning mode	FIX
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H

System - Miscellaneous

Coil Combine Mode	Adaptive Combine
Save uncombined	Off
Matrix Optimization	Off
AutoAlign	---
Coil Select Mode	On - Coil Memory

System - Adjustments

B0 Shim mode	Standard
B1 Shim mode	TrueForm
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto

System - Adjust Volume

! Position	R5.9 A10.9 H6.9 mm
! Orientation	T > C-13.7 > S-3.7
! Rotation	-5.00 deg
! A >> P	210 mm
! R >> L	210 mm
! F >> H	99 mm
Reset	Off

System - pTx Volumes

B1 Shim mode	TrueForm
--------------	----------

System - Tx/Rx

Frequency 1H	123.250023 MHz
Correction factor	1
Gain	High
Img. Scale Cor.	1.000
Reset	Off
? Ref. amplitude 1H	0.000 V

Physio - Signal1

1st Signal/Mode	None
TR	50.00 ms
Concatenations	20

Inline - Common

Subtract	Off
Measurements	325
StdDev	Off
Save original images	On

Inline - MIP

MIP-Sag	Off
MIP-Cor	Off
MIP-Tra	Off
MIP-Time	Off
Save original images	On

Inline - Composing

Inline Composing	Off
Distortion Corr.	Off

Sequence - Part 1

Introduction	Off
Dimension	2D
Contrasts	1
Multi-slice mode	Sequential
Bandwidth	3970 Hz/Px

Sequence - Part 2

Gradient mode	Fast
RF spoiling	Off

Sequence - Nuclei

TX/RX Nucleus	1H
TX/RX delta frequency	0 Hz
TX Nucleus	None
TX delta frequency	0 Hz
Coil elements	HC1-7

Sequence - Special

Reference Scan	Off
Shift	3 FOV/*
Segment	2 #
crusher bfatSAT	0.0 ms
crusher afatSAT	0.4 ms
slewrate	22.00 T/m/s

Sequence - Assistant

Mode	Off
------	-----

\\Research Funded\RESEARCH\FHLIN\leegfmri_emoclips_s015_protocol\ep2d_bold_moco_for_coregistration

TA: 0:36 PM: FIX Voxel size: 3.3×3.3×4.0 mmPAT: Off Rel. SNR: 1.00 : epfid

Properties

Prio recon	Off
Load images to viewer	On
Inline movie	Off
Auto store images	On
Load images to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Auto close inline display	Off
Start measurement without further preparation	Off
Wait for user to start	On
Start measurements	Single measurement

Routine

Slice group	1
Slices	20
Dist. factor	25 %
Position	R5.9 A10.9 H6.9 mm
Orientation	T > C-13.7 > S-3.7
Phase enc. dir.	A >> P
AutoAlign	---
Phase oversampling	0 %
FoV read	210 mm
FoV phase	100.0 %
Slice thickness	4.0 mm
TR	2000 ms
TE	27.0 ms
Averages	1
Concatenations	1
Filter	None
Coil elements	HC1-7

Contrast - Common

TR	2000 ms
TE	27.0 ms
MTC	Off
Flip angle	60 deg
Fat suppr.	Fat sat.

Contrast - Dynamic

Averages	1
Averaging mode	Long term
Reconstruction	Magnitude
Measurements	15
Delay in TR	0 ms
Multiple series	Off

Resolution - Common

FoV read	210 mm
FoV phase	100.0 %
Slice thickness	4.0 mm
Base resolution	64
Phase resolution	100 %
Phase partial Fourier	Off
Interpolation	Off

Resolution - iPAT

Accel. mode	None
-------------	------

Resolution - Filter Image

Distortion Corr.	Off
Prescan Normalize	Off

Resolution - Filter Rawdata

Raw filter	Off
Elliptical filter	Off
Hamming	Off

Geometry - Common

Slice group	1
Slices	20
Dist. factor	25 %
Position	R5.9 A10.9 H6.9 mm
Orientation	T > C-13.7 > S-3.7
Phase enc. dir.	A >> P
FoV read	210 mm
FoV phase	100.0 %
Slice thickness	4.0 mm
TR	2000 ms
Multi-slice mode	Interleaved
Series	Interleaved
Concatenations	1

Geometry - AutoAlign

Slice group	1
Position	R5.9 A10.9 H6.9 mm
Orientation	T > C-13.7 > S-3.7
Phase enc. dir.	A >> P
AutoAlign	---
Initial Position	R5.9 A10.9 H6.9
R	5.9 mm
A	10.9 mm
H	6.9 mm
Initial Rotation	-5.00 deg
Initial Orientation	T > C
T > C	-13.7
> S	-3.7

Geometry - Saturation

Fat suppr.	Fat sat.
Special sat.	None

Geometry - Tim Planning Suite

Set-n-Go Protocol	Off
Table position	H
Table position	0 mm
Inline Composing	Off

System - Miscellaneous

Positioning mode	FIX
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Sum of Squares
Matrix Optimization	Performance
AutoAlign	---

System - Miscellaneous

Coil Select Mode	Off - AutoCoilSelect
------------------	----------------------

System - Adjustments

B0 Shim mode	Standard
B1 Shim mode	TrueForm
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto

System - Adjust Volume

! Position	R5.9 A10.9 H6.9 mm
! Orientation	T > C-13.7 > S-3.7
! Rotation	-5.00 deg
! A >> P	210 mm
! R >> L	210 mm
! F >> H	99 mm
Reset	Off

System - pTx Volumes

B1 Shim mode	TrueForm
Excitation	Standard

System - Tx/Rx

Frequency 1H	123.250023 MHz
Correction factor	1
Gain	High
Img. Scale Cor.	1.000
Reset	Off
? Ref. amplitude 1H	0.000 V

Physio - Signal1

1st Signal/Mode	None
TR	2000 ms
Concatenations	1

BOLD

GLM Statistics	On
Dynamic t-maps	Off
Ignore meas. at start	0
Ignore after transition	0
Model transition states	On
Temp. highpass filter	On
Threshold	4.00
Paradigm size	15
Meas[1]	Baseline
Meas[2]	Baseline
Meas[3]	Baseline
Meas[4]	Baseline
Meas[5]	Baseline
Meas[6]	Baseline
Meas[7]	Baseline
Meas[8]	Baseline
Meas[9]	Baseline
Meas[10]	Baseline
Meas[11]	Baseline
Meas[12]	Baseline
Meas[13]	Baseline
Meas[14]	Baseline
Meas[15]	Active
Motion correction	On
Spatial filter	Off
Measurements	15

BOLD

Delay in TR	0 ms
Multiple series	Off

Sequence - Part 1

Introduction	On
Multi-slice mode	Interleaved
Free echo spacing	Off
Echo spacing	0.49 ms
Bandwidth	2442 Hz/Px

Sequence - Part 2

EPI factor	64
RF pulse type	Normal
Gradient mode	Performance
Excitation	Standard

Sequence - pTX Pulses

\\Research Funded\RESEARCH\FHLIN\eegefMRI_emoclips_s015_protocol\AAHead_Scout

TA: 0:14 PM: REF Voxel size: 1.6×1.6×1.6 mmPAT: 3 Rel. SNR: 1.00 : fl

Properties

Prio recon	On
Load images to viewer	On
Inline movie	Off
Auto store images	On
Load images to stamp segments	Off
Load images to graphic segments	On
Auto open inline display	Off
Auto close inline display	Off
Start measurement without further preparation	Off
Wait for user to start	On
Start measurements	Single measurement

Routine

Slab group	1
Slabs	1
Dist. factor	20 %
Position	L0.0 A10.0 H0.0 mm
Orientation	Sagittal
Phase enc. dir.	A >> P
Phase oversampling	0 %
Slice oversampling	0.0 %
Slices per slab	128
FoV read	260 mm
FoV phase	100.0 %
Slice thickness	1.6 mm
TR	3.15 ms
TE	1.37 ms
Averages	1
Concatenations	1
Filter	Prescan Normalize
Coil elements	HC1-7;NC1,2

Contrast - Common

TR	3.15 ms
TE	1.37 ms
Flip angle	8 deg

Contrast - Dynamic

Averages	1
Averaging mode	Short term
Reconstruction	Magnitude
Measurements	1

Resolution - Common

FoV read	260 mm
FoV phase	100.0 %
Slice thickness	1.6 mm
Base resolution	160
Phase resolution	100 %
Slice resolution	69 %
Phase partial Fourier	6/8
Slice partial Fourier	6/8
Trajectory	Cartesian

Resolution - iPAT

PAT mode	GRAPPA
Accel. factor PE	3
Ref. lines PE	24
Accel. factor 3D	1

Resolution - iPAT

Reference scan mode	Integrated
---------------------	------------

Resolution - Filter Image

Image Filter	Off
Distortion Corr.	Off
Prescan Normalize	On
Unfiltered images	Off
Normalize	Off
B1 filter	Off

Resolution - Filter Rawdata

Raw filter	Off
Elliptical filter	Off

Geometry - Common

Slab group	1
Slabs	1
Dist. factor	20 %
Position	L0.0 A10.0 H0.0 mm
Orientation	Sagittal
Phase enc. dir.	A >> P
Slice oversampling	0.0 %
Slices per slab	128
FoV read	260 mm
FoV phase	100.0 %
Slice thickness	1.6 mm
TR	3.15 ms
Multi-slice mode	Sequential
Series	Ascending
Concatenations	1

Geometry - AutoAlign

Slab group	1
Position	L0.0 A10.0 H0.0 mm
Orientation	Sagittal
Phase enc. dir.	A >> P
Initial Position	Isocenter
L	0.0 mm
P	0.0 mm
H	0.0 mm
Initial Rotation	0.00 deg
Initial Orientation	Transversal

Geometry - Tim Planning Suite

Set-n-Go Protocol	Off
Table position	H
Table position	0 mm
Inline Composing	Off

System - Miscellaneous

Positioning mode	REF
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Adaptive Combine
Save uncombined	Off
Matrix Optimization	Off

System - Miscellaneous

Coil Select Mode	Off - AutoCoilSelect
------------------	----------------------

System - Adjustments

B0 Shim mode	Tune up
B1 Shim mode	TrueForm
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto

System - Adjust Volume

Position	Isocenter
Orientation	Transversal
Rotation	0.00 deg
A >> P	263 mm
R >> L	350 mm
F >> H	350 mm
Reset	Off

System - pTx Volumes

B1 Shim mode	TrueForm
Excitation	Non-sel.

System - Tx/Rx

Frequency 1H	123.250023 MHz
Correction factor	1
Gain	Low
Img. Scale Cor.	1.000
Reset	Off
? Ref. amplitude 1H	0.000 V

Physio - PACE

Resp. control	Off
Concatenations	1

Inline - Common

Flip angle	8 deg
Measurements	1
Time to center	6.2 s

Inline - Inline

Subtract	Off
Measurements	1
StdDev	Off
Save original images	On

Inline - MIP

MIP-Sag	Off
MIP-Cor	Off
MIP-Tra	Off
MIP-Time	Off
Save original images	On

Inline - Composing

Inline Composing	Off
Distortion Corr.	Off

Inline - MapIt

Save original images	On
MapIt	None
Flip angle	8 deg
Measurements	1

Inline - MapIt

Contrasts	1
TR	3.15 ms
TE	1.37 ms

Sequence - Part 1

Introduction	On
Dimension	3D
Asymmetric echo	Weak
Contrasts	1
Multi-slice mode	Sequential
Bandwidth	540 Hz/Px

Sequence - Part 2

RF pulse type	Fast
Gradient mode	Normal
Excitation	Non-sel.
RF spoiling	On

Sequence - Assistant

Mode	Off
------	-----

\\Research Funded\RESEARCH\FHLIN\eeegfmri_emoclips_s015_protocol\MBSIREPI_ref

TA: 5.0 s PM: FIX Voxel size: 5.0×5.0×4.0 mmRel. SNR: 1.00 : fl

Properties

Prio recon	Off
Load images to viewer	On
Inline movie	Off
Auto store images	On
Load images to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Auto close inline display	Off
Start measurement without further preparation	Off
Wait for user to start	On
Start measurements	Single measurement

Routine

Slice group	1
Slices	20
Dist. factor	25 %
Position	R5.9 A10.9 H6.9 mm
Orientation	T > C-13.7 > S-3.7
Phase enc. dir.	A >> P
AutoAlign	---
Phase oversampling	0 %
FoV read	210 mm
FoV phase	100.0 %
Slice thickness	4.0 mm
TR	50.00 ms
TE	27.60 ms
Averages	1
Concatenations	20
Filter	None
Coil elements	HC1-7

Contrast - Common

TR	50.00 ms
TE	27.60 ms
TD	0 ms
MTC	Off
Flip angle	30 deg
Fat suppr.	Fat sat.
Water suppr.	None

Contrast - Dynamic

Averages	1
Averaging mode	Short term
Reconstruction	Magnitude
Measurements	10
Pause after meas. 1	0.0 s
Pause after meas. 2	0.0 s
Pause after meas. 3	0.0 s
Pause after meas. 4	0.0 s
Pause after meas. 5	0.0 s
Pause after meas. 6	0.0 s
Pause after meas. 7	0.0 s
Pause after meas. 8	0.0 s
Pause after meas. 9	0.0 s
Multiple series	Off

Resolution - Common

FoV read	210 mm
FoV phase	100.0 %

Resolution - Common

Slice thickness	4.0 mm
Base resolution	42
Phase resolution	100 %
Phase partial Fourier	Off
Interpolation	Off

Resolution - Filter Image

Image Filter	Off
Distortion Corr.	Off
Prescan Normalize	Off
Normalize	Off
B1 filter	Off

Resolution - Filter Rawdata

Raw filter	Off
Elliptical filter	Off

Geometry - Common

Slice group	1
Slices	20
Dist. factor	25 %
Position	R5.9 A10.9 H6.9 mm
Orientation	T > C-13.7 > S-3.7
Phase enc. dir.	A >> P
FoV read	210 mm
FoV phase	100.0 %
Slice thickness	4.0 mm
TR	50.00 ms
Multi-slice mode	Sequential
Series	Ascending
Concatenations	20

Geometry - AutoAlign

Slice group	1
Position	R5.9 A10.9 H6.9 mm
Orientation	T > C-13.7 > S-3.7
Phase enc. dir.	A >> P
AutoAlign	---
Initial Position	R5.9 A10.9 H6.9
R	5.9 mm
A	10.9 mm
H	6.9 mm
Initial Rotation	-5.00 deg
Initial Orientation	T > C
T > C	-13.7
> S	-3.7

Geometry - Saturation

Fat suppr.	Fat sat.
Water suppr.	None
Special sat.	None

Geometry - Tim Planning Suite

Set-n-Go Protocol	Off
Table position	H
Table position	0 mm
Inline Composing	Off

System - Miscellaneous

Positioning mode	FIX
------------------	-----

System - Miscellaneous

Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Adaptive Combine
Save uncombined	Off
Matrix Optimization	Off
AutoAlign	---
Coil Select Mode	On - Coil Memory

System - Adjustments

B0 Shim mode	Standard
B1 Shim mode	TrueForm
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto

System - Adjust Volume

! Position	R5.9 A10.9 H6.9 mm
! Orientation	T > C-13.7 > S-3.7
! Rotation	-5.00 deg
! A >> P	210 mm
! R >> L	210 mm
! F >> H	99 mm
Reset	Off

System - pTx Volumes

B1 Shim mode	TrueForm
--------------	----------

System - Tx/Rx

Frequency 1H	123.250023 MHz
Correction factor	1
Gain	High
Img. Scale Cor.	1.000
Reset	Off
? Ref. amplitude 1H	0.000 V

Physio - Signal1

1st Signal/Mode	None
TR	50.00 ms
Concatenations	20

Inline - Common

Subtract	Off
Measurements	10
StdDev	Off
Save original images	On

Inline - MIP

MIP-Sag	Off
MIP-Cor	Off
MIP-Tra	Off
MIP-Time	Off
Save original images	On

Inline - Composing

Inline Composing	Off
Distortion Corr.	Off

Sequence - Part 1

Introduction	Off
Dimension	2D
Contrasts	1
Multi-slice mode	Sequential
Bandwidth	3970 Hz/Px

Sequence - Part 2

Gradient mode	Fast
RF spoiling	Off

Sequence - Nuclei

TX/RX Nucleus	1H
TX/RX delta frequency	0 Hz
TX Nucleus	None
TX delta frequency	0 Hz
Coil elements	HC1-7

Sequence - Special

Reference Scan	On
Shift	3 FOV/*
Segment	2 #
crusher bfatSAT	0.0 ms
crusher afatSAT	0.4 ms
slewrate	22.00 T/m/s

Sequence - Assistant

Mode	Off
------	-----

\\Research Funded\RESEARCH\FHLIN\eeegfmri_emoclips_s015_protocol\emoclips3_MBSIREPI_acc_1.4s

TA: 8:00 PM: FIX Voxel size: 5.0×5.0×4.0 mmRel. SNR: 1.00 : fl

Properties

Prio recon	Off
Load images to viewer	On
Inline movie	Off
Auto store images	On
Load images to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Auto close inline display	Off
Start measurement without further preparation	Off
Wait for user to start	On
Start measurements	Single measurement

Routine

Slice group	1
Slices	20
Dist. factor	25 %
Position	R5.9 A10.9 H6.9 mm
Orientation	T > C-13.7 > S-3.7
Phase enc. dir.	A >> P
AutoAlign	---
Phase oversampling	0 %
FoV read	210 mm
FoV phase	100.0 %
Slice thickness	4.0 mm
TR	50.00 ms
TE	27.60 ms
Averages	1
Concatenations	20
Filter	None
Coil elements	HC1-7

Contrast - Common

TR	50.00 ms
TE	27.60 ms
TD	0 ms
MTC	Off
Flip angle	30 deg
Fat suppr.	Fat sat.
Water suppr.	None

Contrast - Dynamic

Averages	1
Averaging mode	Short term
Reconstruction	Magnitude
Measurements	321
Pause after meas.	1.4 s
Multiple series	Off

Resolution - Common

FoV read	210 mm
FoV phase	100.0 %
Slice thickness	4.0 mm
Base resolution	42
Phase resolution	100 %
Phase partial Fourier	Off
Interpolation	Off

Resolution - Filter Image

Image Filter	Off
Distortion Corr.	Off
Prescan Normalize	Off
Normalize	Off
B1 filter	Off

Resolution - Filter Rawdata

Raw filter	Off
Elliptical filter	Off

Geometry - Common

Slice group	1
Slices	20
Dist. factor	25 %
Position	R5.9 A10.9 H6.9 mm
Orientation	T > C-13.7 > S-3.7
Phase enc. dir.	A >> P
FoV read	210 mm
FoV phase	100.0 %
Slice thickness	4.0 mm
TR	50.00 ms
Multi-slice mode	Sequential
Series	Ascending
Concatenations	20

Geometry - AutoAlign

Slice group	1
Position	R5.9 A10.9 H6.9 mm
Orientation	T > C-13.7 > S-3.7
Phase enc. dir.	A >> P
AutoAlign	---
Initial Position	R5.9 A10.9 H6.9
R	5.9 mm
A	10.9 mm
H	6.9 mm
Initial Rotation	-5.00 deg
Initial Orientation	T > C
T > C	-13.7
> S	-3.7

Geometry - Saturation

Fat suppr.	Fat sat.
Water suppr.	None
Special sat.	None

Geometry - Tim Planning Suite

Set-n-Go Protocol	Off
Table position	H
Table position	0 mm
Inline Composing	Off

System - Miscellaneous

Positioning mode	FIX
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H

System - Miscellaneous

Coil Combine Mode	Adaptive Combine
Save uncombined	Off
Matrix Optimization	Off
AutoAlign	---
Coil Select Mode	On - Coil Memory

System - Adjustments

B0 Shim mode	Standard
B1 Shim mode	TrueForm
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto

System - Adjust Volume

! Position	R5.9 A10.9 H6.9 mm
! Orientation	T > C-13.7 > S-3.7
! Rotation	-5.00 deg
! A >> P	210 mm
! R >> L	210 mm
! F >> H	99 mm
Reset	Off

System - pTx Volumes

B1 Shim mode	TrueForm
--------------	----------

System - Tx/Rx

Frequency 1H	123.250023 MHz
Correction factor	1
Gain	High
Img. Scale Cor.	1.000
Reset	Off
? Ref. amplitude 1H	0.000 V

Physio - Signal1

1st Signal/Mode	None
TR	50.00 ms
Concatenations	20

Inline - Common

Subtract	Off
Measurements	321
StdDev	Off
Save original images	On

Inline - MIP

MIP-Sag	Off
MIP-Cor	Off
MIP-Tra	Off
MIP-Time	Off
Save original images	On

Inline - Composing

Inline Composing	Off
Distortion Corr.	Off

Sequence - Part 1

Introduction	Off
Dimension	2D
Contrasts	1
Multi-slice mode	Sequential
Bandwidth	3970 Hz/Px

Sequence - Part 2

Gradient mode	Fast
RF spoiling	Off

Sequence - Nuclei

TX/RX Nucleus	1H
TX/RX delta frequency	0 Hz
TX Nucleus	None
TX delta frequency	0 Hz
Coil elements	HC1-7

Sequence - Special

Reference Scan	Off
Shift	3 FOV/*
Segment	2 #
crusher bfatSAT	0.0 ms
crusher afatSAT	0.4 ms
slewrate	22.00 T/m/s

Sequence - Assistant

Mode	Off
------	-----

\\Research Funded\RESEARCH\FHLIN\leegfmri_emoclips_s015_protocol\ep2d_bold_moco_for_coregistration

TA: 0:36 PM: FIX Voxel size: 3.3×3.3×4.0 mmPAT: Off Rel. SNR: 1.00 : epfid

Properties

Prio recon	Off
Load images to viewer	On
Inline movie	Off
Auto store images	On
Load images to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Auto close inline display	Off
Start measurement without further preparation	Off
Wait for user to start	On
Start measurements	Single measurement

Routine

Slice group	1
Slices	20
Dist. factor	25 %
Position	R5.9 A10.9 H6.9 mm
Orientation	T > C-13.7 > S-3.7
Phase enc. dir.	A >> P
AutoAlign	---
Phase oversampling	0 %
FoV read	210 mm
FoV phase	100.0 %
Slice thickness	4.0 mm
TR	2000 ms
TE	27.0 ms
Averages	1
Concatenations	1
Filter	None
Coil elements	HC1-7

Contrast - Common

TR	2000 ms
TE	27.0 ms
MTC	Off
Flip angle	60 deg
Fat suppr.	Fat sat.

Contrast - Dynamic

Averages	1
Averaging mode	Long term
Reconstruction	Magnitude
Measurements	15
Delay in TR	0 ms
Multiple series	Off

Resolution - Common

FoV read	210 mm
FoV phase	100.0 %
Slice thickness	4.0 mm
Base resolution	64
Phase resolution	100 %
Phase partial Fourier	Off
Interpolation	Off

Resolution - iPAT

Accel. mode	None
-------------	------

Resolution - Filter Image

Distortion Corr.	Off
Prescan Normalize	Off

Resolution - Filter Rawdata

Raw filter	Off
Elliptical filter	Off
Hamming	Off

Geometry - Common

Slice group	1
Slices	20
Dist. factor	25 %
Position	R5.9 A10.9 H6.9 mm
Orientation	T > C-13.7 > S-3.7
Phase enc. dir.	A >> P
FoV read	210 mm
FoV phase	100.0 %
Slice thickness	4.0 mm
TR	2000 ms
Multi-slice mode	Interleaved
Series	Interleaved
Concatenations	1

Geometry - AutoAlign

Slice group	1
Position	R5.9 A10.9 H6.9 mm
Orientation	T > C-13.7 > S-3.7
Phase enc. dir.	A >> P
AutoAlign	---
Initial Position	R5.9 A10.9 H6.9
R	5.9 mm
A	10.9 mm
H	6.9 mm
Initial Rotation	-5.00 deg
Initial Orientation	T > C
T > C	-13.7
> S	-3.7

Geometry - Saturation

Fat suppr.	Fat sat.
Special sat.	None

Geometry - Tim Planning Suite

Set-n-Go Protocol	Off
Table position	H
Table position	0 mm
Inline Composing	Off

System - Miscellaneous

Positioning mode	FIX
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Sum of Squares
Matrix Optimization	Performance
AutoAlign	---

System - Miscellaneous

Coil Select Mode	Off - AutoCoilSelect
------------------	----------------------

System - Adjustments

B0 Shim mode	Standard
B1 Shim mode	TrueForm
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto

System - Adjust Volume

! Position	R5.9 A10.9 H6.9 mm
! Orientation	T > C-13.7 > S-3.7
! Rotation	-5.00 deg
! A >> P	210 mm
! R >> L	210 mm
! F >> H	99 mm
Reset	Off

System - pTx Volumes

B1 Shim mode	TrueForm
Excitation	Standard

System - Tx/Rx

Frequency 1H	123.250023 MHz
Correction factor	1
Gain	High
Img. Scale Cor.	1.000
Reset	Off
? Ref. amplitude 1H	0.000 V

Physio - Signal1

1st Signal/Mode	None
TR	2000 ms
Concatenations	1

BOLD

GLM Statistics	On
Dynamic t-maps	Off
Ignore meas. at start	0
Ignore after transition	0
Model transition states	On
Temp. highpass filter	On
Threshold	4.00
Paradigm size	15
Meas[1]	Baseline
Meas[2]	Baseline
Meas[3]	Baseline
Meas[4]	Baseline
Meas[5]	Baseline
Meas[6]	Baseline
Meas[7]	Baseline
Meas[8]	Baseline
Meas[9]	Baseline
Meas[10]	Baseline
Meas[11]	Baseline
Meas[12]	Baseline
Meas[13]	Baseline
Meas[14]	Baseline
Meas[15]	Active
Motion correction	On
Spatial filter	Off
Measurements	15

BOLD

Delay in TR	0 ms
Multiple series	Off

Sequence - Part 1

Introduction	On
Multi-slice mode	Interleaved
Free echo spacing	Off
Echo spacing	0.49 ms
Bandwidth	2442 Hz/Px

Sequence - Part 2

EPI factor	64
RF pulse type	Normal
Gradient mode	Performance
Excitation	Standard

Sequence - pTX Pulses

\\Research Funded\RESEARCH\FHLIN\eeegfmri_emoclips_s015_protocol\t1_mprage_sag_p2 iso1.0

TA: 8:02 PM: REF Voxel size: 1.0×1.0×1.0 mmPAT: 2 Rel. SNR: 1.00 : tfl

Properties

Prio recon	Off
Load images to viewer	On
Inline movie	Off
Auto store images	On
Load images to stamp segments	On
Load images to graphic segments	On
Auto open inline display	Off
Auto close inline display	Off
Start measurement without further preparation	Off
Wait for user to start	On
Start measurements	Single measurement

Routine

Slab group	1
Slabs	1
Dist. factor	50 %
Position	R4.2 A16.5 F18.2 mm
Orientation	S > C-3.1 > T1.1
Phase enc. dir.	A >> P
AutoAlign	Head > Brain
Phase oversampling	100 %
Slice oversampling	45.5 %
Slices per slab	176
FoV read	256 mm
FoV phase	100.0 %
Slice thickness	1.00 mm
TR	1800.0 ms
TE	2.21 ms
Averages	1
Concatenations	1
Filter	Prescan Normalize
Coil elements	HC1-7

Contrast - Common

TR	1800.0 ms
TE	2.21 ms
Magn. preparation	Non-sel. IR
T1	900 ms
Flip angle	10 deg
Fat suppr.	None
Water suppr.	None

Contrast - Dynamic

Averages	1
Averaging mode	Long term
Reconstruction	Magnitude
Measurements	1
Multiple series	Each measurement

Resolution - Common

FoV read	256 mm
FoV phase	100.0 %
Slice thickness	1.00 mm
Base resolution	256
Phase resolution	100 %
Slice resolution	100 %
Phase partial Fourier	Off
Slice partial Fourier	Off
Interpolation	Off

Resolution - iPAT

PAT mode	GRAPPA
Accel. factor PE	2
Ref. lines PE	24
Accel. factor 3D	1
Reference scan mode	Integrated

Resolution - Filter Image

Image Filter	Off
Distortion Corr.	Off
Prescan Normalize	On
Unfiltered images	Off
Normalize	Off
B1 filter	Off

Resolution - Filter Rawdata

Raw filter	Off
Elliptical filter	Off

Geometry - Common

Slab group	1
Slabs	1
Dist. factor	50 %
Position	R4.2 A16.5 F18.2 mm
Orientation	S > C-3.1 > T1.1
Phase enc. dir.	A >> P
Slice oversampling	45.5 %
Slices per slab	176
FoV read	256 mm
FoV phase	100.0 %
Slice thickness	1.00 mm
TR	1800.0 ms
Multi-slice mode	Single shot
Series	Ascending
Concatenations	1

Geometry - AutoAlign

Slab group	1
Position	R4.2 A16.5 F18.2 mm
Orientation	S > C-3.1 > T1.1
Phase enc. dir.	A >> P
AutoAlign	Head > Brain
Initial Position	L1.6 A1.3 F28.2
L	1.6 mm
A	1.3 mm
F	28.2 mm
Initial Rotation	-0.07 deg
Initial Orientation	S > T
S > T	-2.0
> C	1.2

Geometry - Navigator**Geometry - Tim Planning Suite**

Set-n-Go Protocol	Off
Table position	H
Table position	0 mm
Inline Composing	Off

System - Miscellaneous

Positioning mode	REF
------------------	-----

System - Miscellaneous

Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Adaptive Combine
Save uncombined	Off
Matrix Optimization	Off
AutoAlign	Head > Brain
Coil Select Mode	Off - All

System - Adjustments

B0 Shim mode	Standard
B1 Shim mode	TrueForm
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto

System - Adjust Volume

Position	R4.2 A16.5 F18.2 mm
Orientation	S > C-3.1 > T1.1
Rotation	3.17 deg
A >> P	256 mm
F >> H	256 mm
R >> L	176 mm
Reset	Off

System - pTx Volumes

B1 Shim mode	TrueForm
Excitation	Non-sel.

System - Tx/Rx

Frequency 1H	123.250023 MHz
Correction factor	1
Gain	Low
Img. Scale Cor.	1.000
Reset	Off
? Ref. amplitude 1H	0.000 V

Physio - Signal1

1st Signal/Mode	None
TR	1800.0 ms
Concatenations	1

Physio - Cardiac

Magn. preparation	Non-sel. IR
Ti	900 ms
Fat suppr.	None
Dark blood	Off
FoV read	256 mm
FoV phase	100.0 %
Phase resolution	100 %

Physio - PACE

Resp. control	Off
Concatenations	1

Inline - Common

Subtract	Off
Measurements	1

Inline - Common

StdDev	Off
Save original images	On

Inline - MIP

MIP-Sag	Off
MIP-Cor	Off
MIP-Tra	Off
MIP-Time	Off
Save original images	On

Inline - Composing

Inline Composing	Off
Distortion Corr.	Off

Inline - MapIt

Save original images	On
MapIt	None
Flip angle	10 deg
Measurements	1
TR	1800.0 ms
TE	2.21 ms

Sequence - Part 1

Introduction	On
Dimension	3D
Elliptical scanning	Off
Reordering	Linear
Asymmetric echo	Allowed
Flow comp.	No
Multi-slice mode	Single shot
Echo spacing	6.8 ms
Bandwidth	210 Hz/Px

Sequence - Part 2

RF pulse type	Normal
Gradient mode	Performance
Excitation	Non-sel.
RF spoiling	On
Incr. Gradient spoiling	On
Turbo factor	256

Sequence - Assistant

Mode	Off
------	-----