

## Introduction

This manual is edited to help researchers interested in using magnetic resonance imaging (MRI) data in MEG source modeling. We will specifically use the following software packages.

### I. FreeSurfer

An open source software suite for processing and analyzing (human) brain MRI images. For more information, please check

<https://surfer.nmr.mgh.harvard.edu/>

fMRI data pre-processing basically involves the following steps:

- Unpack DICOM files into the file format allowing us for subsequent processing
- Motion correction
- Slice-timing correction
- Spatial smoothing

## Environment setup

- **Login**

1. Open a SSH session
2. Login to server

```
%ssh 140.119.165.24 -l username1 -Y
```

```
(%ssh -o TCPKeepAlive=no -o ServerAliveInterval=15  
140.119.165.24 -l username -Y)
```

3. Goto your working directory<sup>2</sup>

```
> cd  
/space/maki5/1/users/fhlin/ini_vm_nccu/10102012_ChangCI/epi_data
```

4. Prepare environment

```
> source /space/maki/1/pubsw/bme-dev-env-dev.csh
```

---

<sup>1</sup> Here I use “user\_name” as the example.

<sup>2</sup> Here the working directory is  
/space/maki5/1/users/fhlin/ini\_vm\_nccu/10102012\_ChangCI/epi\_data.  
You should change that accordingly in your own analysis.

## Unpack DICOM files

- **Scan EPI sessions**

The first step is to check how many EPI “sessions” in the data stored in the folder VM\_EPI.

```
> unpacksdcmdir -src VM_EPI -targ . -scanonly ./info
```

By the end of screen output, you may find the following:

```
Done scanning Wed Feb 19 16:47:01 CST 2014
-----
  4          ep2d_bold_MGH  ok   64  64  30 120
10102012_CHANGCI.MR.FHLIN_FMRI.0004.0001.2012.10.10.19.21.1
1.750000.5743858.IMA
```

It indicates there is only ONE (1) session in this data. And the EPI session has session number 4.

- **Perform unpacking**

Now ONE (1) file is needed to specify which run(s) is going to be unpacked. Here I created a file called unpack.rule. The content of the file is:

```
> cat unpack.rule
4    bold nii  f.nii
```

It specifies that run 4 should be unpacked inside the folder BOLD with nii<sup>3</sup> file format. The output file should be bold/004/f.nii.

Now perform the unpacking with this configuration rule.

```
> unpacksdcmdir -src VM_EPI -targ . -cfg ./unpack.rule
```

If things run smoothly, you would see the last few lines of your screen showing this:

```
StartTime: Wed Feb 19 16:51:22 CST 2014
EndTime:   Wed Feb 19 16:51:35 CST 2014
unpacksdcmdir Done
```

Double-check if files are there as expected:

---

<sup>3</sup> Nii is the file name extension for the NiftI file format, a common file format storing neuroimaging data. Please see <http://nifti.nimh.nih.gov/nifti-1/> for details.

```
> ls -al bold/004
total 78424
drwxrwxr-x 2 fhlin fhlin      4096 Apr 24  2014 .
drwxrwxr-x 3 fhlin fhlin      4096 Feb 19 16:56 ..
-rw-rw-r-- 1 fhlin fhlin      9600 Apr 24  2014 flf
-rw-rw-r-- 1 fhlin fhlin 29491552 Apr 24  2014 f.nii
-rw-rw-r-- 1 fhlin fhlin    69458 Apr 24  2014 f.nii-
infodump.dat
```

It specifies that run 4 should be unpacked inside the folder BOLD with nii file format. The output file should be bold/004/f.nii.

## Pre-processing steps

Two files are needed to specify all folders containing unpacked EPI data. Here I created `sessid` and `sessdir` two texts files for this purpose:

```
>pwd
/autofs/space/maki5_001/users/fhlin/ini_vm_nccu/10102012_Ch
angCI/epi_data

> ls
rolle7:/space/maki5/1/users/fhlin/ini_vm_nccu/10102012_Ch
angCI/epi_data> ls -al
total 120
drwxr-xr-x  5 fhlin fhlin  4096 Apr 24 17:52 .
drwxr-xr-x 11 fhlin fhlin  4096 Apr 24 17:50 ..
drwxrwxr-x  3 fhlin fhlin  4096 Feb 19 16:56 bold
-rw-rw-r--  1 fhlin fhlin 18000 Apr 24 17:52
dicomdir.sumfile
-rw-rw-r--  1 fhlin fhlin   126 Feb 19 16:47 info
drwxrwxr-x  2 fhlin fhlin  4096 Feb 19 21:49 log
-rw-rw-r--  1 fhlin fhlin    56 Feb 19 16:33 sessdir
-rw-rw-r--  1 fhlin fhlin     9 Feb 19 16:33 sessid
-rw-rw-r--  1 fhlin fhlin     0 Feb 19 21:47 stc-session
-rw-rw-r--  1 fhlin fhlin 50544 Apr 24 17:52 unpack.log
-rw-rw-r--  1 fhlin fhlin    17 Feb 19 16:49 unpack.rule
drwxr-xr-x  2 fhlin fhlin 16384 Dec  4 01:52 VM_EPI

>cat sessid
epi_data

>cat sessdir
/space/maki5/1/users/fhlin/ini_vm_nccu/10102012_ChangCI
```

- **Make analysis/pre-processing template**

```
> mktemplate-sess -sf sessid -df sessdir
```

- **Motion correction**

All files inside each folder within 'bold' will be motion corrected.

```
> mc-sess -sf sessid -df sessdir -per-run
```

Check the output files:

```
>ls -l bold/004
total 79904
-rw-rw-r-- 1 fhlin fhlin      9600 Apr 24 17:52 flf
-rw-rw-r-- 1 fhlin fhlin    20210 Apr 24 2014
fmcpr.mat.aff12.1D
-rw-rw-r-- 1 fhlin fhlin     11040 Apr 24 17:54 fmcpr.mcdat
-rw-rw-r-- 1 fhlin fhlin 52054603 Apr 24 2014 fmcpr.nii.gz
-rw-rw-r-- 1 fhlin fhlin      5362 Apr 24 2014
fmcpr.nii.gz.mclog
-rw-rw-r-- 1 fhlin fhlin 29491552 Apr 24 17:52 f.nii
-rw-rw-r-- 1 fhlin fhlin     69458 Apr 24 17:52 f.nii-
infodump.dat
-rw-rw-r-- 1 fhlin fhlin       470 Apr 24 2014
mcdat2extreg.log
-rw-rw-r-- 1 fhlin fhlin      6480 Apr 24 2014 mcprextreg
-rw-rw-r-- 1 fhlin fhlin        21 Apr 24 2014 template.log
-rw-rw-r-- 1 fhlin fhlin    119787 Apr 24 2014
template.nii.gz
drwxrwxr-x 3 fhlin fhlin      4096 Apr 24 2014
tmpdir.mcdat2extreg.31315
drwxrwxr-x 3 fhlin fhlin      4096 Apr 24 2014
tmpdir.mcdat2extreg.31428
```

The file `fmcpr.nii.gz` contains the motion-corrected EPI data in zipped nii format.

- **Slice timing correction**

Slice timing correction is going to apply to files with file stem 'fmcpr'. The output file after slice timing correction has the file stem 'fmcprstc'.

```
> stc-sess -sf sessid -df sessdir -i fmcpr -o fmcprstc
```

Check the output files:

```
> ls -l bold/004
total 131328
-rw-rw-r-- 1 fhlin fhlin      9600 Apr 24 17:52 flf
-rw-rw-r-- 1 fhlin fhlin    20210 Apr 24 17:56
fmcpr.mat.aff12.1D
-rw-rw-r-- 1 fhlin fhlin     11040 Apr 24 17:54 fmcpr.mcdat
-rw-rw-r-- 1 fhlin fhlin 52054603 Apr 24 17:56 fmcpr.nii.gz
-rw-rw-r-- 1 fhlin fhlin      5362 Apr 24 17:56
fmcpr.nii.gz.mclog
```

```

-rw-rw-r-- 1 fhlin fhlin 52653539 Apr 24 2014
fmcprstc.nii.gz
-rw-rw-r-- 1 fhlin fhlin      3211 Apr 24 2014
fmcprstc.nii.gz.log
-rw-rw-r-- 1 fhlin fhlin 29491552 Apr 24 17:52 f.nii
-rw-rw-r-- 1 fhlin fhlin  69458 Apr 24 17:52 f.nii-
infodump.dat
-rw-rw-r-- 1 fhlin fhlin      470 Apr 24 17:56
mcdat2extreg.log
-rw-rw-r-- 1 fhlin fhlin      6480 Apr 24 17:56 mcprextreg
-rw-rw-r-- 1 fhlin fhlin       21 Apr 24 17:55 template.log
-rw-rw-r-- 1 fhlin fhlin  119787 Apr 24 17:55
template.nii.gz
drwxrwxr-x 3 fhlin fhlin      4096 Apr 24 17:55
tmpdir.mcdat2extreg.31315
drwxrwxr-x 3 fhlin fhlin      4096 Apr 24 17:55
tmpdir.mcdat2extreg.31428

```

The file `fmcprstc.nii.gz` contains the motion-corrected AND slice timing corrected EPI data in zipped nii format.

- **Spatial smoothing**

Spatial smoothing is going to apply to files with file stem 'fmcprstc' using a Gaussian smooth kernel with full-width-half-maximum (FWHM) of 10 mm. The output file after slice timing correction has the file stem 'fmcprstcs'.

```
> spatialsMOOTH-session -sf sessid -df sessdir -i fmcprstc -o
fmcprstcs -fwhm 10 -outfmt nii -no-mask
```

Check the output files:

```

> ls -l bold/004
total 188932
-rw-rw-r-- 1 fhlin fhlin      9600 Apr 24 17:52 flf
-rw-rw-r-- 1 fhlin fhlin     20210 Apr 24 17:56
fmcpr.mat.aff12.1D
-rw-rw-r-- 1 fhlin fhlin      11040 Apr 24 17:54 fmcpr.mcdat
-rw-rw-r-- 1 fhlin fhlin 52054603 Apr 24 17:56 fmcpr.nii.gz
-rw-rw-r-- 1 fhlin fhlin      5362 Apr 24 17:56
fmcpr.nii.gz.mclog
-rw-rw-r-- 1 fhlin fhlin 52653539 Apr 24 2014
fmcprstc.nii.gz
-rw-rw-r-- 1 fhlin fhlin      3211 Apr 24 2014
fmcprstc.nii.gz.log
-rw-rw-r-- 1 fhlin fhlin 58982752 Apr 24 2014
fmcprstcs.nii

```

```
-rw-rw-r-- 1 fhlin fhlin 29491552 Apr 24 17:52 f.nii
-rw-rw-r-- 1 fhlin fhlin 69458 Apr 24 17:52 f.nii-
infodump.dat
-rw-rw-r-- 1 fhlin fhlin 470 Apr 24 17:56
mcdat2extreg.log
-rw-rw-r-- 1 fhlin fhlin 6480 Apr 24 17:56 mcprextreg
-rw-rw-r-- 1 fhlin fhlin 21 Apr 24 17:55 template.log
-rw-rw-r-- 1 fhlin fhlin 119787 Apr 24 17:55
template.nii.gz
drwxrwxr-x 3 fhlin fhlin 4096 Apr 24 17:55
tmpdir.mcdat2extreg.31315
drwxrwxr-x 3 fhlin fhlin 4096 Apr 24 17:55
tmpdir.mcdat2extreg.31428
```

The file `fmcpdstc.nii` contains the motion-corrected, slice timing corrected, AND spatially smoothed EPI data in nii format.